

## HOW TO MEASURE FOR TributeNight™ Leg & Lower Torso Garments

①

Mark the length reference points from where you would like the garment to start and end at the determined positions directly on your patient's skin using a water-based marker.

②

Measure the circumferences in centimeters at those same marked reference points and record the measurements on the left side of the diagram on your Order Form. *(See details below.)*

③

Measure the length at each marked reference point, and record the measurements on the right side of the diagram on your Order Form. *(See details below.)*

*We recommend patients be in a relaxed, supine (laying down) position.*

### Answers to Frequently Asked Questions

In extreme presentations, measure the lobules, breast and/or pannus in a supported position.

When skin folds are present, measure the smallest and largest circumferences.

Contact our Custom Design Center to send photos or illustrations to further explain patient presentations.

With patients who fluctuate frequently, please feel free to add a circumference range.

### Circumferences

Left side of diagram

- A<sup>c</sup>** Measure foot circumference at determined **A** position (distal end of garment on foot).
- B<sup>c</sup>** Measure circumference at determined **B** position (largest ankle circumference).
- Y<sup>c</sup>** Measure circumference at determined **Y** position (heel/zero position) by taking a diagonal measurement over anterior foot back to posterior heel.
- C<sup>c</sup>** Measure circumference of calf at determined **C** position (largest calf circumference).
- D<sup>c</sup>** Measure circumference just below patella at determined **D** position (below patella).
- E<sup>c</sup>** Measure circumference at determined **E** position (mid-patella).
- F<sup>c</sup>** Measure circumference of upper leg at determined **F** position (largest thigh circumference).
- G<sup>c</sup>** Measure circumference straight across at determined **G** position (groin).
- H<sup>c</sup>** Measure circumference across hips at determined **H** position (Iliac Crest).
- I<sup>c</sup>** Measure circumference at determined **I** position (natural waistline).
- J<sup>c</sup>** Measure circumference at determined **J** position (under breast).

### Lengths

Right side of diagram

- A<sup>L</sup>** Measure foot length from **Y** (heel/zero position) to **A** (distal end of garment on foot).
- B<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **B** (largest ankle circumference).
- C<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **C** (largest calf circumference).
- D<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **D** (below patella).
- E<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **E** (mid-patella).
- F<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **F** (largest thigh circumference).
- LG<sup>L</sup>** Optional: Measure lateral leg length from **Y** (heel/zero position) to **G** (groin). *A pencil can be placed in groin as a visual guide.*
- MG<sup>L</sup>** Measure medial leg length from medial **Y** (zero position) to medial **G** length (groin).
- H<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **H** (Iliac Crest).
- I<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **I** (natural waistline).
- J<sup>L</sup>** Measure lateral leg length from **Y** (heel/zero position) to **J** (under chest).
- AS<sup>L</sup>** Measure **AS** (anterior straddle length) from mid-crotch to desired anterior proximal end of garment.
- PS<sup>L</sup>** Measure **PS** (posterior straddle length) from mid-crotch to desired posterior proximal end of garment.