



# FOOT COMPRESSION WEAR



STANDARD & CUSTOM FOOT COMPRESSION WEAR



## COMPREBOOT™ - LITE

- Ultra-thin profile - easily fits into most footwear.
- Highly adjustable design, may be trimmed if needed.
- Breathable stretch fabric conforms well to ankle area.
- Provides mild, therapeutic compression.

SMALL - LARGE		XL/XXL	
REGULAR	LONG	REGULAR	LONG
1101-FL-REG	1101-FL-LNG	1102-FL-REG	1102-FL-LNG



## COMPREBOOT™

- Features neoprene material with hook/loop enclosers.
- Open heel & toe.
- Provides moderate, therapeutic compression.

	SMALL		MED/LARGE		XL/XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
BLACK	1101-F-REG	1101-F-LNG	1102-F-REG	1102-F-LNG	1103-F-REG	1103-F-LNG
BEIGE	1111-F-REG	1111-F-LNG	1112-F-REG	1112-F-LNG	1113-F-REG	1113-F-LNG



## COMPREBOOT™ PLUS

- Chip foam lining, gripper patches for traction and D-ring straps.
- Straps are removeable for custom positioning or replacement.

	SMALL		MED/LARGE		XL/XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
	1101-FP-REG	1101-FP-LNG	1102-FP-REG	1102-FP-LNG	1103-FP-REG	1103-FP-LNG




## MEDABOOT™ - CUSTOM

- Ideal for non-standard foot sizes.
- Added dorsum padding.
- Features a breathable, moisture wicking liner.
- Photos required.

## MEASURING INSTRUCTIONS:

- 1: Record **length** measurement **K** (length from heel to 1st metatarsal head, or desired length) on the corresponding line. Select the boot size that matches according to **row K** on the sizing chart.
- 1b: (**MedaBoot only**) Record **length** measurements **M** (top of foot/3rd Matatarsal Head to Ankle Bend) on the corresponding lines.
- 2: Record **circumference** measurements **I, J, and L** on the corresponding lines.
- 3: Check to make certain that foot measurements **I, J, and L** do **not exceed** the adjustment range of the boot size selected in **Step 1** (measurement **K** - sizing chart **row K**).
- 4: If any of the above foot measurements (**I, J, and L**) exceed the adjustment range of the boot size selected in **Step 1** (**K**) above, select the **next larger size** where all three measurements do not exceed the adjustment range of the boot.



	SMALL		MED/LARGE		XL/XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
<b>I</b>	28 - 36	28 - 36	39 max	39 max	44 max	44 max
<b>J</b>	22 - 26	22 - 26	30 max	30 max	33 max	33 max
<b>K</b>	14 - 18	19 - 23	18 - 20	23 - 25	20 - 22	25 - 27
	25 - 30	25 - 30	36 max	36 max	43 max	43 max