

STANDARD HAND COMPRESSION WEAR



MEDAHAND™

- Breathable/moisture wicking/anti-odor and soil-release fabric.
- Open cell foam chips provide gentle massaging effect.
- Conforming design ensures proper fit and gradient compression.
- OverSleeve available for additional compression (pg. 65).

	SMALL	MEDIUM	LARGE	XLARGE
RIGHT	1201-H-R	1202-H-R	1203-H-R	1204-H-R
LEFT	1201-H-L	1202-H-L	1203-H-L	1204-H-L



MEDAGLOVE™

- Open cell foam chips provide gentle massaging effect.
- Conforming design with finger separation ensures proper fit.
- Effective as a compression garment, or bandage liner.
- Gradient compression from finger tips to wrist.
- OverSleeve available for additional compression (pg. 65).


	SMALL	MEDIUM	LARGE	XLARGE
RIGHT	1201-G-R	1202-G-R	1203-G-R	1204-G-R
LEFT	1201-G-L	1202-G-L	1203-G-L	1204-G-L



DORSAL POCKET GLOVE

- Dorsal pocket to insert foam or other compression product.
- Full finger compression with doffing loops for easy use.
- Palm pad to fill in palm space.
- Ultra-thin profile for full hand function during use.

	SMALL	MEDIUM	LARGE	XLARGE
RIGHT	1201-DG-R	1202-DG-R	1203-DG-R	1204-DG-R
LEFT	1201-DG-L	1202-DG-L	1203-DG-L	1204-DG-L



HAND SIZING CHART

Use this chart to measure for hand/glove products. Item numbers listed under the garment description.

	SMALL	MEDIUM	LARGE	XLARGE
PALM (E)	15 - 19	19 - 22	22 - 26	26 - 29
WRIST (A)	14 - 20	16 - 22	18 - 24	20 - 26