



CIRCAID REDUCTION KIT LEG GARMENT MEASUREMENT FORM

Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

upper leg kit: Check the box next to the corresponding length and width.

***cG cannot be more than 30cm larger than cE2.**

| | Length | | | Width | | Quantity |
|-------|----------------------|-------------------------|--------------------|-------------------------|---------------------------|----------|
| | Short | Standard | Long | Regular | Wide | |
| Left | ∠ KC-Groin < 37cm | ∠ KC-Groin 37 - 41cm | ∠ KC-Groin > 41 | cG < 80cm cE2 < 65cm | cG < 125cm cE2 < 110cm | |
| Right | | | | | | |

Each upper leg kit includes:

- 1 reduction upper leg component
- 1 reduction shelf strap
- 6 Velcro[®] stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use.

***cE2 cannot be more than 30cm larger than cC.**

| | Length | Width | Quantity |
|-------|--------|--|----------|
| Left | 30cm | cE2 < 110cm cE < 102cm cC < 80cm | |
| Right | | | |

Each knee kit includes:

- 1 reduction knee component
- 2 Velcro[®] knee spines
- 1 paper measuring tape
- 1 direction for use
- 1 reduction shelf strap

lower leg kit: Check the box next to the corresponding length and width.

***cC cannot be more than 30cm larger than cB.**

| | Length | | | Width | | Quantity |
|-------|---------------------|------------------------|---------------------|------------------------|------------------------|----------|
| | Short | Standard | Long | Regular | Wide | |
| Left | ∠ mall-KC < 32cm | ∠ mall-KC 32 - 36cm | ∠ mall-KC > 36cm | cC < 65cm cB < 55cm | cC < 80cm cB < 65cm | |
| Right | | | | | | |

kit includes:

- 1 reduction kit undersleeve leg
- 1 reduction shelf strap
- 6 Velcro[®] stays
- 1 paper measuring tape
- 1 direction for use

foot options: Check the box next to the corresponding length and width.

*** foot options are not included with the kit and must be ordered separately.**

| | compressionanklet | | pacband™ | Single band AFW [®] | Quantity |
|-------|-------------------|-------------|-----------|------------------------------|-------------|
| | Standard | Large | X-Large | | |
| Left | cH 19-36 cm | cH 22-41 cm | cH < 44cm | cH > 41cm | |
| Right | | | | | cB 19-27 cm |

additional ancillary items:

| reduction kit undersleeve leg Quantity (Pairs) | reduction kit shelf strap Quantity (Each) |
|---|--|
| | |

