



Account Number: 120167

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

# CUSTOM LEGGING SIZING CHART


CIRCUMFERENCE MEASUREMENTS		LANDMARK MEASURING POINTS
LEFT	RIGHT	
		<b>E</b> circumference at center of patella*
		<b>B</b> ankle circumference around malleoli
		5cm
		10cm
		15cm
		20cm
		25cm
		30cm
		35cm
		40cm
		45cm
		50cm
		55cm
		60cm
		65cm
		70cm
		75cm
		80cm
		85cm

- Record all measurements in centimeters. Take circumference measurements starting at the malleolus and every 5cm above. For Lower Leggings end circumference measurements at knee crease. All other garments continue taking circumference measurements up the leg, including the patella, to the groin (required for proper knee-and-above positioning and placement).
- Complete length and foot measurements. Provide contour and linear lengths when measuring an uncommon limb shape.
- Take measurements in shaded areas for standard sized leggings and ankle-foot wraps. NOTE: 1in = 2.54cm

## JUXTA-FIT™ STANDARD LEGGING SIZES

	Small	Medium	Medium Full Calf	Large	Large Full Calf	XL	XL Full Calf	XXL
5cm	15-23cm	20-28cm	20-28cm	25-33cm	25-33cm	30-38cm	30-38cm	33-41cm
10cm								
15cm	20-28cm	25-33cm	33-40cm	30-38cm	38-46cm	33-41cm	41-51cm	43-53cm
20cm								
25cm	30-38cm	36-43cm	43-50cm	41-48cm	48-56cm	41-51cm	51-61cm	53-64cm
30cm								
35cm								
40cm								
45cm								
50cm								
55cm								
60cm								
65cm								
70cm								
75cm								
80cm								
85cm								

LENGTH MEASUREMENTS		
LEFT	RIGHT	
		<b>B-D</b> medial length from prominent point of malleolus to knee crease
		<36cm Short >36cm Long
		<b>B-E</b> lateral length from prominent point of malleolus to center of patella*
		<b>B-G</b> medial length from prominent point of malleolus to groin, far enough below pubic area so that garment will not cause discomfort*
		<b>B-K</b> lateral length from prominent point of malleolus to gluteal fold*




**MEASUREMENT STARTS FROM MALLEOLUS, NOT FROM FLOOR**  
 NOTE: Measurements B-D & B-G are where garment will typically begin and end. If Shelves, Folds, or Trouble Areas are Present, Please Draw Shape of Patient's Leg Darker on Sample or on Separate Sheet.  
 \*Required for ALL except Lower Legging

PRODUCT SELECTION	
<input type="checkbox"/>	28cm Juxta-Fit™ Premium Standard Lower Legging (Short) <sup>5</sup>
<input type="checkbox"/>	36cm Juxta-Fit™ Premium Standard Lower Legging (Long) <sup>5</sup>
<input type="checkbox"/>	28cm Juxta-Fit™ Essentials Standard Lower Legging (Short) <sup>5</sup>
<input type="checkbox"/>	36cm Juxta-Fit™ Essentials Standard Lower Legging (Long) <sup>5</sup>
<input type="checkbox"/>	Juxta-Fit™ Custom Lower Legging (only) <sup>4</sup>
<input type="checkbox"/>	Juxta-Fit™ Custom Knee Piece (only)
<input type="checkbox"/>	Juxta-Fit™ Custom Upper Legging (only)
<input type="checkbox"/>	Juxta-Fit™ Custom Lower Legging w/Knee <sup>4</sup>
<input type="checkbox"/>	Juxta-Fit™ Custom Upper Legging w/Knee
<input type="checkbox"/>	Juxta-Fit™ Whole Legging (please choose CoverUp color below)
<input type="checkbox"/>	<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
<input type="checkbox"/>	Classic Flex™ Custom Lower Legging
<input type="checkbox"/>	Graduate™ Whole Leg (please choose CoverUp color below)
<input type="checkbox"/>	<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
<input type="checkbox"/>	Graduate™ BK
<input type="checkbox"/>	Graduate™ Boot Only

GRADUATE™ COLOR OPTIONS				
Exterior Foam Liner	<input type="checkbox"/> BLACK	<input type="checkbox"/> BEIGE	<input type="checkbox"/> RED	<input type="checkbox"/> TURQUOISE
Interior Foam Liner	N/A	<input type="checkbox"/> BEIGE	<input type="checkbox"/> RED	<input type="checkbox"/> TURQUOISE

SELECT FOOT OPTION <sup>4</sup>	
<input type="checkbox"/>	Comfort Compression Ankle™ (max H circumference 34cm)
<input type="checkbox"/>	Comfort EZ Single-Band Ankle Foot Wrap™ (3")
<input type="checkbox"/>	Comfort EZ Single-Band Ankle Foot Wrap™ (4")
<input type="checkbox"/>	Standard Juxta-Fit™ Ankle-Foot Wrap (max H circumference 34cm)
<input type="checkbox"/>	Standard Juxta-Fit™ Interlocking Ankle-Foot Wrap (max H circ. 34cm)
<input type="checkbox"/>	Custom Juxta-Fit™ Ankle Foot Wrap
<input type="checkbox"/>	<input type="checkbox"/> ATTACHED <input type="checkbox"/> SEPARATE
<input type="checkbox"/>	Juxta-Lite™ Ankle-Foot Wrap (max H circumference 31cm)

FOOT MEASUREMENTS		LANDMARK MEASURING POINTS	JUXTA-FIT™ & JUXTA-LITE™ STANDARD ANKLE FOOT WRAP (AFW) SIZES		
LEFT	RIGHT		S	M	L
		<b>B</b> ankle circumference around malleoli	< 36cm	< 36cm	< 41cm
		<b>H</b> arch circumference	19-24cm	24-29cm	29-34cm
		<b>A-Y</b> length from base of great toe to where foot meets ankle	> 9cm	> 9cm	> 9cm
		<b>Z</b> length from base of toes to back of heel (For Interlocking AFW)	13-25cm	13-25cm	13-25cm
		<b>Z</b> length for Standard Juxta-Fit™ & Juxta-Lite™ AFW <sup>4</sup>	> 16cm	> 16cm	> 16cm
		<b>A-B</b> length from floor to center of medial malleoli			
		<b>A</b> circumference at base of toes			
		<b>Y</b> circumference around heel to where foot meets ankle			
		<b>P<sub>2</sub></b> length from malleoli to malleoli, pass under heel			



<sup>4</sup> EZ Single-Band Standard Ankle-Foot Wrap™ included by default  
<sup>5</sup> Comfort Compression Anklets™ included by default