



Patient Last Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_  
 Date: \_\_\_\_\_

Account Number: 120167  
 Patient First Name: \_\_\_\_\_  
 Fitter First Name: \_\_\_\_\_  
 (example PT/OT/PTA)



## Non-Custom Lower Leg Garment Measurement Form & Sizing Chart

### SIZING CHART

Patient should be standing when taking measurements. Limbs should be at their smallest; limb size should be maintained from time of measurements to assure a proper fit of garment. When on borderline between sizes, select a smaller size if significant edema is present and reduction is anticipated; if edema is minimal, select larger size. After measuring the smallest dimension of the ankle and the largest dimension of the calf, if any of the measurements do not fall within one of these size ranges, call (800) 247-2243. Length measurement is for Juxta-Fit™ only, actual legging lengths are 28cm (Short) and 36cm (Long); Ready-Fit™ and Classic-Flex™ have one adjustable length. **Note 1in = 2.54cm**

MEASUREMENTS		
	LEFT	RIGHT
<b>B</b>		
<b>B<sub>1</sub></b>		
<b>C</b>		
<b>H</b>		
<b>B-D</b>		

CLASSIC-FLEX™ & READY-FIT™			
SIZE	CIRCUMFERENCE MEASUREMENT		
	B	B <sub>1</sub>	C
SMALL	15-23cm	20-28cm	30-38cm
MEDIUM	20-28cm	25-33cm	36-43cm
LARGE	25-33cm	30-38cm	41-48cm
XL	30-38cm	33-41cm	41-51cm
XXL	33-41cm	43-53cm	53-64cm

JUXTA-FIT™ PREMIUM & ESSENTIALS			
SIZE	CIRCUMFERENCE MEASUREMENT		
	B	B <sub>1</sub>	C
SMALL	15-23cm	20-28cm	30-38cm
MEDIUM	20-28cm	25-33cm	36-43cm
MEDIUM FULLCALF	20-28cm	33-40cm	43-50cm
LARGE	25-33cm	30-38cm	41-48cm
LARGE FULLCALF	25-33cm	38-46cm	48-56cm
XL	30-38cm	33-41cm	41-51cm
XL FULLCALF	30-38cm	41-51cm	43-53cm
XXL*	33-41cm	43-53cm	53-64cm

AFW MEASUREMENT	
SIZE	H
SMALL	19-24cm
MEDIUM	24-29cm
LARGE	29-34cm


LENGTH MEASUREMENT	
SIZE	B-D
SHORT	under 36cm
LONG	over 36cm

\* If measurements fall within XXL size range, please consider ordering a custom Juxta-Fit™ for optimal fit.  
 † EZ Single-Band Ankle-Foot Wrap™ included by default  
 ‡ Juxta-Lite™ Ankle-Foot Wrap is available only in sizes Small and Medium  
 § Comfort Compression Anklets™ included by default

LEGGING OPTIONS	SIZE:	QTY:
28CM JUXTA-FIT™ PREMIUM † STANDARD LOWER LEGGING (SHORT)		
36CM JUXTA-FIT™ PREMIUM † STANDARD LOWER LEGGING (LONG)		
28CM JUXTA-FIT™ ESSENTIALS † STANDARD LOWER LEGGING (SHORT)		
36CM JUXTA-FIT™ ESSENTIALS † STANDARD LOWER LEGGING (LONG)		
READY-FIT™ STANDARD LOWER LEGGING ‡		
CLASSIC-FLEX™ STANDARD LOWER LEGGING ‡		

FOOT OPTIONS	SIZE:	QTY:
Comfort Compression Anklets™ (max H circumference 34cm)		
Comfort EZ Single-Band Ankle Foot Wrap™ (3")		
Comfort EZ Single-Band Ankle Foot Wrap™ (4")		
Standard Juxta-Fit™ Ankle-Foot Wrap (max H circumference 34cm)		
Standard Juxta-Fit™ Interlocking Ankle-Foot Wrap (max H circ. 34cm)		
Custom Juxta-Fit™ Interlocking Ankle-Foot Wrap		
Juxta-Lite™ Ankle-Foot Wrap (max H circumference 31cm) §		

ANCILLARIES	SIZE:	QTY:
COMFORT COVERUP™ FOR LOWER LEG (BLACK) (Small: 34cm - 50cm)		
COMFORT COVERUP™ FOR LOWER LEG (BLACK) (Large: 51cm - 80cm)		
COMFORT COVERUP™ FOR LOWER LEG (BEIGE) (Small: 34cm - 50cm)		
COMFORT COVERUP™ FOR LOWER LEG (BEIGE) (Large: 51cm - 80cm)		
COMFORT CAPRI™ (For a sizing chart visit www.circaid.com)		
INSTRUCTIONAL DVD		



PRICES SUBJECT TO CHANGE WITHOUT NOTICE.