




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



FARROWRAP OTS CALF LENGTH LEG GARMENT MEASUREMENT FORM & SIZING CHART

Note: All measurements in cm



All FarrowWrap Off-The-Shelf legpieces are sized by taking three simple measurements: widest calf circumference (C), ankle circumference (B), and floor to knee crease length (A-D).

All FarrowWrap Off-The-Shelf footpieces are sized by taking two simple measurements: mid-foot circumference (A¹) and straight distance from back of heel to base of fifth metatarsal head (X).

Legpiece					
	X-Small	Small	Medium	Large	X-Large
(C) Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
(B) Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
(A-D) Regular	33-36 cm	35-38 cm	37-40 cm	39-42 cm	39-42 cm
(A-D) Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm

Footpiece					
	X-Small	Small	Medium	Large	X-Large
(A¹) Mid-Foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
(X) Regular	16-17.5 cm	17.5-19 cm	19-20.5cm	20.5-22 cm	22.5-24 cm
(X) Long	18-19.5 cm	19.5-21 cm	21-22.5 cm	22.5-24 cm	24.5-26 cm



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



FARROWRAP OTS CALF LENGTH LEG GARMENT MEASUREMENT FORM & SIZING CHART

	Classic	LITE	STRONG	BASIC	ThinGrip*
Legpiece Regular XS					*
Legpiece Regular S					*
Legpiece Regular M					*
Legpiece Regular L					*
Legpiece Regular XL					*
 					
Legpiece Tall XS					*
Legpiece Tall S					*
Legpiece Tall M					*
Legpiece Tall L					*
Legpiece Tall XL					*
 					
Footpiece Regular XS					
Footpiece Regular S					
Footpiece Regular M					
Footpiece Regular L					
Footpiece Regular XL					
 					
Footpiece Long XS					
Footpiece Long S					
Footpiece Long M					
Footpiece Long L					
Footpiece Long XL					

* ThinGrip Off-The-Shelf available in Footpiece only