




**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

## FarrowWrap® 4000 Knee High Order Form

Note: All measurements in cm



Please place a numeric value in the white box to indicate items needed.	FarrowWrap® 4000 Legpiece Size Chart							
	XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty
(C) Widest Calf	25 – 38 cm		32 – 48 cm		37 – 56 cm		43 – 67 cm	
(B) Ankle	17 – 23 cm		19 – 28 cm		25 – 37 cm		29 – 43 cm	
(B-D) Regular	30 – 33 cm		30 – 33 cm		30 – 33 cm		30 – 33 cm	
(B-D) Tall	34 – 37 cm		34 – 37 cm		34 – 37 cm		34 – 37 cm	

Color Options	
Tan	<input type="checkbox"/>
Black	<input type="checkbox"/>

*Default color option is Tan. A matching FarrowHybrid® will be included with the FarrowWrap®4000*

### FarrowWrap® 4000 Legpiece

Obtain the following 3 measurements:

- Widest calf circumference (at **C**)
- Least ankle circumference (at **B**) just above malleoli.
- Length of leg following posterior contour from the least ankle to 2 finger widths below knee crease, **B-D**

### Note

Each FarrowWrapN 4000 legpiece comes with one (1) color matching FarrowHybrid® Foot Compression Liner. **FarrowHybrid® AD Foot Compression** is part compression stocking (over the foot and ankle) and part sock (ankle to knee) and may be used in place of our footpieces with our legpiece. Please be sure to indicate the color of the FarrowWrap® 4000 desired. If no color is selected, the default color will be tan.

**THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY**