



C

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____



Lower Extremity Order Form

IMPORTANT - All measurements must be recorded in centimeters. Style must be selected. Please write clear!
 Elvarex Soft is available in Class 2 (23-32 mmHg), Beige only in Styles AD (Knee High), AG (Thigh High) and AT (Waist High) all with Slant Open or Slant Closed Toe, Silicone band on top only.

Quantity / Class	Left	Right	Panty
Elvarex 1 (18-21 mmHg*)			
Elvarex 2 (23-32 mmHg*)			
Elvarex Soft 2 (23-32 mmHg*)			
Elvarex 3 (36-46 mmHg*)			
Elvarex-forte 3F (36-46 mmHg*)			
Elvarex 4 (59-70 mmHg*)			
Elvarex-super 4S (70-90 mmHg*)			

Color

- BEIGE - Elvarex Soft Available in Beige only
- BLACK Available Class 1-4 Only

Circum. (c)	Length (l)		
cT	K2-T	lT	
cH	K1-T	lH	
Circumference (c)	Length (l): Taken from each landmark to floor		
	Left	Right	Right
cG		lG	
cF		lF	
cE		lE	
cD		lD	
cC		lC	
cB1		lB1	
cB		lB	
cY		lA (medial)	
cA		lA (lateral)	

* Design Pressure

CAUTION : This product contains natural rubber latex which may cause allergic reactions.



Special Styles

- AD Knee High
- AF Mid Thigh
- AG Thigh High
- AGTL Chap Left
- AGTR Chap Right
- AGT Chap Pair (Both Legs)
- AT Pantyhose
- AT Pantyhose, 1 1/2 leg
- GT Biker Shorts
- FG Leg Extension

All lengths are linear
 Leg /Total lengths - from each landmark to floor

Options

- Silicone Dotted Band 2.5cm
 - Top
 - Inside
- Silicone Dotted Band 5cm
 - Top
 - Inside
- Zipper B to D only
 - Inside (Medial)
 - Outside (Lateral)
 - Seam
- Zipper E to G only
 - Inside (Medial)
 - Outside (Lateral)
 - Seam
- Fly for Men Open Pubis
- Adjustable Waistband
- Straight Open Toe _____ cm (Lateral only)
- Straight Closed Toe _____ cm (Total length)
- Slant Closed Toe _____ cm (Total length)
- T-Heel (Available in Class 2-3F Only)



Comments:
