




Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_



Quantity/Class	CCL 1 18-21 mmHg*	CCL 2 23-32 mmHg*	CCL 3 34-46 mmHg*	CCL 3F 34-46 mmHg*	CCL 4 49-70 mmHg*	CCL 4S 60-90 mmHg*
Left						
Right						
Pressure panty Elvarex (Body Bandage)						

# ELVAREX® Lower Extremity Order Form

Circum. (c)	Length (l)	Length (l)
cT	K2-T	lT
cH	K1-T	lH



### Basic styles

- AD Knee
- AF Mid-Thigh
- AG Thigh
- AG-T (AG with chapstyle)
  - Piece
  - Pair
- AG-HT 1½ Leg Pantyhose
- AT Pantyhose

### Options

- Profile/Ankle pad
- T-Heel (Class 2-3 forte only)

### Comments

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### Variations

- B'G
- BG
- FG (leg extension)
- GT Biker Shorts
- B1G-T Capri

### Color


- Beige
- Black
- Cocoa
- Navy
- Grey

### Special options

- SoftFit™ (Only Available CCL 1-3)
- Silicone dotted band 2.5cm:
  - Top
  - Inside
  - Inside ¾
- Silicone dotted band 5cm:
  - Top
  - Inside
  - Inside ¾
- Vertical Silicone Strips AG\*\*\*:
  - Front
  - Back
  - Both
- Zipper B to D only:
  - Inside (medial)
  - Outside (Lateral)
- Zipper E to G only:
  - Inside (medial)
  - Outside (Lateral)
- Knee Comfort Zone (Not available in ccl 1)
- Pocket
  - Instep (Not Available with Profile)
  - Back of Knee\*\*\*
  - All four sides closed
- Fly for Men
- Open Pubis
- Adjustable Waistband
- Top Comfort Zone

Circumference (c)		Length (l): Taken from each landmark to floor	
Left	Right	Left	Right
cG		lG	
cF		lF	
cE		lE	
cD		lD	
cC		lC	
cB1		lB1	
cB		lB	
cY		lA (medial)	
cA		lA (lateral)	

- Straight Open Toe Lateral Length \_\_\_\_\_ cm
- Slant Open Toe Medial Length \_\_\_\_\_ cm
- Slant Closed Toe Medial Length \_\_\_\_\_ cm
- Straight Closed Toe Total Foot Length \_\_\_\_\_ cm
- Lateral Length \_\_\_\_\_ cm
- Total Foot Length \_\_\_\_\_ cm



\* Design Pressure  
 \*\* Only available in knee high.  
 \*\*\* Not available with knee comfort.  
 CAUTION: This product contains natural rubber latex which may cause allergic reactions.

**THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY**