



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

JOBST NON-CUSTOM CIRCULAR KNIT ARM SLEEVE MEASUREMENT FORM AND SIZING CHART



Measure in centimeters

Left Right

Length Measurement from axilla to wrist (c-g) = _____

LONG LENGTH IS > 43 CMS REGULAR LENGTH IS < 43 CMS*

The length measurement is taken from **Wrist** along the inside of the arm to the inseam, **Axilla**, of the arm; measure with patient standing or sitting with arm down by the side.

Measurement Points

Circumference **Axilla** is taken around the upper arm at the axilla. For ease in determining this point, place a book up into the arm pit area. Point 'g' on the arm will be even with the top of the book.

Circumference **Wrist** is the point of greatest compression and therefore the most important point. This point is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point. **Be sure not to measure too tight.**

Circumference Wrist = _____ with Silicone* with no Silicone

Circumference 'Mid-Lower' = _____

Circumference 'Mid-Upper' = _____

Circumference Axilla = _____

Material	Class I 15-20 mmHg	Class II 20-30 mmHg	Class II 30-40 mmHg	Colors Available
Bella Strong	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>
Bella Lite Arm Sleeve	<input type="checkbox"/>	<input type="checkbox"/>		Beige <input type="checkbox"/> Brown <input type="checkbox"/> Black <input type="checkbox"/> * Silicone only available in Beige * Brown and Black are only available in 20-30 mmHg *Brown and Black are only available in LONG length

Sizing Chart - Bella Lite Armsleeve

Sizing Chart - Bella Strong Armsleeve

Jobst Bella Lite Arm Sleeve

	WRIST	MID-LOWER	MID-UPPER
S	5 ³ / ₈ " - 7 ⁵ / ₈ " (14-20 cm)	7" - 9 ¹ / ₈ " (18-23 cm)	8 ¹ / ₂ " - 15 ¹ / ₂ " (22-40 cm)
M	6 ¹ / ₈ " - 8 ¹ / ₂ " (16-22 cm)	9 ¹ / ₄ " - 11 ¹ / ₄ " (23-29 cm)	10 ⁵ / ₈ " - 17 ³ / ₄ " (27-45 cm)
L	6 ⁷ / ₈ " - 9 ¹ / ₄ " (18-24 cm)	11 ³ / ₈ " - 13 ¹ / ₂ " (29-34 cm)	12 ³ / ₄ " - 19 ⁷ / ₈ " (33-51 cm)

Lengths: Regular: up to 17" (43 cm); Long: over 17" (43 cm)



Size	Wrist		Mid-Lower		Top Axilla ¹	
	Min.	Max.	Min.	Max.	Min.	Max.
1	5 ¹ / ₂ " (14 cm)	7" (18 cm)	6 ³ / ₄ " (17 cm)	8 ¹ / ₄ " (21 cm)	8 ¹ / ₂ " (22 cm)	12 ¹ / ₄ " (31 cm)
2	5 ¹ / ₂ " (14 cm)	7" (18 cm)	6 ³ / ₄ " (17 cm)	8 ¹ / ₄ " (21 cm)	12 ¹ / ₄ " (31 cm)	16 ¹ / ₄ " (41 cm)
3	5 ¹ / ₂ " (14 cm)	7" (18 cm)	8 ¹ / ₄ " (21 cm)	9 ³ / ₄ " (25 cm)	9 ³ / ₄ " (25 cm)	13 ¹ / ₄ " (34 cm)
4	6 ¹ / ₄ " (16 cm)	7 ¹ / ₂ " (19 cm)	8 ¹ / ₄ " (21 cm)	9 ³ / ₄ " (25 cm)	11" (28 cm)	15" (38 cm)
5	6 ³ / ₄ " (17 cm)	7 ³ / ₄ " (20 cm)	8 ¹ / ₄ " (21 cm)	9 ³ / ₄ " (25 cm)	13 ¹ / ₄ " (34 cm)	17" (43 cm)
6	5 ¹ / ₂ " (14 cm)	7" (18 cm)	9 ³ / ₄ " (25 cm)	11 ¹ / ₂ " (29 cm)	11 ¹ / ₂ " (29 cm)	15" (38 cm)
7	6 ³ / ₄ " (17 cm)	7 ³ / ₄ " (20 cm)	9 ³ / ₄ " (25 cm)	11 ¹ / ₂ " (29 cm)	13" (33 cm)	17" (43 cm)
8	7" (18 cm)	8 ¹ / ₄ " (21 cm)	9 ³ / ₄ " (25 cm)	11 ¹ / ₂ " (29 cm)	15" (38 cm)	19" (48 cm)
9	6 ¹ / ₄ " (16 cm)	7 ¹ / ₂ " (19 cm)	11 ¹ / ₂ " (29 cm)	13" (33 cm)	12 ¹ / ₂ " (32 cm)	16 ¹ / ₂ " (42 cm)
10	7 ¹ / ₂ " (19 cm)	8 ¹ / ₂ " (22 cm)	11 ¹ / ₂ " (29 cm)	13" (33 cm)	14 ¹ / ₄ " (36 cm)	18" (46 cm)

¹Axilla measurement to be taken in vertical position. ¹La medida de la axilla debe tomarse en posición vertical.
¹L'axilla doit être mesurée en position verticale.

THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY