



Account Number: 4057807

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_



# Lower Extremity Order Form

**IMPORTANT** - All measurements must be recorded in centimeters. Style must be selected. Please write clearly. Elvarex Soft is available in Class 2 (23-32 mmHg), Beige only in Styles AD (Knee High), AG (Thigh High) and AT (Waist High) all with Slant Open or Slant Closed Toe, Silicone band on top only.

Quantity / Class	Left	Right	Panty
Elvarex 1 (18-21 mmHg*)			
Elvarex 2 (23-32 mmHg*)			
Elvarex Soft 2 (23-32 mmHg*)			
Elvarex 3 (36-46 mmHg*)			
Elvarex-forte 3F (36-46 mmHg*)			
Elvarex 4 (59-70 mmHg*)			
Elvarex-super 4S (70-90 mmHg*)			

### Special Styles

- AD Knee High
- AF Mid Thigh
- AG Thigh High
- AGTL Chap Left
- AGTR Chap Right
- AGT Chap Pair (Both Legs)
- AT Pantyhose
- AT Pantyhose, 1 1/2 leg
- GT Biker Shorts
- FG Leg Extension

All lengths are linear  
 Leg / Total lengths - from each landmark to floor



### Options

- Silicone Dotted Band 2.5cm
  - Top
  - Inside
- Silicone Dotted Band 5cm
  - Top
  - Inside
- Zipper B to D only
  - Inside (Medial)
  - Outside (Lateral)
  - Seam
- Zipper E to G only
  - Inside (Medial)
  - Outside (Lateral)
  - Seam
- Fly for Men  Open Pubis
- Adjustable Waistband
  - Straight Open Toe \_\_\_\_\_ cm (Lateral only)
  - Straight Closed Toe \_\_\_\_\_ cm (Total length)
  - Slant Closed Toe \_\_\_\_\_ cm (Total length)
- T-Heel (Available in Class 2-3F Only)

### Color

- BEIGE - Elvarex Soft - Available in Beige only
- BLACK - Available Class 1-4 Only

Circum. (c)	Length (l)	Length (l)
cT	K2-T	IT
cH	K1-T	IH
Circumference (c)		Length (l): <small>Taken from each landmark to floor</small>
Left	Right	Left Right
cG		IG
cF		IF
cE		IE
cD		ID
cC		IC
cB1		IB1
cB		IB
cY		IA (medial)
cA		IA (lateral)



**TIPS:**  
 Straight open toe - Lateral Length  
 Straight closed toe - Total Foot Length  
 Slant closed toe - Total, Medial & Lateral Length

### Comments:

\* Design Pressure

CAUTION: This product contains natural rubber latex which may cause allergic reactions.