



Account Number: 4057807

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

JOBST NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASUREMENT FORM AND SIZING CHART 

Measure in Centimeters

Left Right

Thigh Length _____

Calf Length _____

A-g REGULAR LENGTH >68 CM
 A-g PETITE LENGTH <68 CM

A-d REGULAR LENGTH IS >38 CM
 A-d PETITE LENGTH IS <38 CM

For thigh length, the length measurement is taken from the floor (at the inside of heel) along the inside of the leg to the inseam/groin, point 'g,' this length should end at the gluteal fold.

For calf length, the length measurement is taken from the floor (at the inside of heel) along the inside of the leg to approximately 3 1/2 centimeters below the crease of the knee, point 'd'

Silicone borders are not usually recommended for calf length stockings, however, they can be necessary for cone-shaped legs.

For calf and thigh length stockings, measure with the patient standing.

Measurement Points

Circumference **point 'g'** is found at the inseam/groin

Circumference **point 'c'** is the largest circumference of the calf

Circumference **point 'b'** is the point of greatest compression and therefore the most important point. This point is the smallest circumference of the ankle found closest to and above the ankle bone

Circumference 'g' = _____


Circumference 'c' = _____

Circumference 'b' = _____

Open Toe Silicone Closed Toe No Silicone


Measurement Point for Pantyhose

Hip circumference = _____



Material	Class 1 15-20 mmHg	Class 2 20-30 mmHg	Class 3 30-40	Colors Available
Relief (calf length) AD	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>
Relief (thigh length) AG	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>
Relief (thigh length) AG		<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>
Relief (pantyhose)			<input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/>

Sizing chart

SIZE	AN	C	THI	
S	7"-8 1/4" (18-21 cm)	11"-15" (28-38 cm)	15 3/4"-24 3/8" (40-62 cm)	28"-46" (71-117 cm)
M	8 7/8"-9 7/8" (21-25 cm)	11 7/8"-16 1/2" (30-42 cm)	18 1/8"-27 1/2" (46-70 cm)	30"-50" (76-127 cm)
L	10"-11 3/8" (25-29 cm)	12 1/2"-18 1/8" (32-46 cm)	21 1/4"-30 3/4" (54-78 cm)	32"-54" (81-137 cm)
XL	11 1/2"-13" (29-33 cm)	13 3/8"-19 5/8" (34-50 cm)	23 5/8"-32" (60-81 cm)	40"-65" (102-165 cm)
LFC	10"-11 7/8" (25-30 cm)	18"-24" (46-61 cm)		
XLFC	12"-14" (30-36 cm)	18"-24" (46-61 cm)		
	Regular: Greater than 15" (38cm). Available in all above sizes. Petite: 15" (38cm) or less. *M & L available in petite sizes.			

Knee High: Measure ankle and calf. Thigh High: Measure ankle, calf and thigh. Waist High and Chaps: Measure ankle, thigh and hip.
 LFC: Large Full Calf XLFC: Extra-Large Full Calf.