



Account Number: 4057807

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

## JOBST NON-CUSTOM HAND GLOVE/GAUNTLET MEASUREMENT FORM AND SIZING CHART



### Measure in centimeters

Left  Right

### Measurement Points:

Circumference **point 'a'** is found across the metacarpus.  
 Circumference **point 'c'** is found distally on the arm where the wrist meets the hand.  
 One should find the smallest circumference at this point.

Circumference 'a' = \_\_\_\_\_

Circumference 'c' = \_\_\_\_\_

Style	
Gauntlet w/ thumb stub only	<input type="checkbox"/>
Glove w/ finger stubs	<input type="checkbox"/>

Material	Class 1 15-20 mmHg	Class 2 20-30 mmHg	Color Available
Bella Strong Gauntlet		<input type="checkbox"/>	Natural <input type="checkbox"/> Honey <input type="checkbox"/>
Ready-to-Wear Gauntlet	<input type="checkbox"/>	<input type="checkbox"/>	Beige <input type="checkbox"/> Brown <input type="checkbox"/> Black <input type="checkbox"/> *Brown and Black only available in 20-30 mmHg

### Sizing chart:

#### Ready-to-Wear Gauntlet

Palm



5<sup>7</sup>/<sub>8</sub>" - 7<sup>3</sup>/<sub>8</sub>" (15cm - 19cm)  
 7<sup>1</sup>/<sub>2</sub>" - 8<sup>5</sup>/<sub>8</sub>" (19cm - 22cm)  
 8<sup>3</sup>/<sub>4</sub>" - 10<sup>1</sup>/<sub>4</sub>" (22cm - 26cm)

Wrist



5<sup>3</sup>/<sub>8</sub>" - 7<sup>5</sup>/<sub>8</sub>" (14cm - 20cm)  
 6<sup>1</sup>/<sub>8</sub>" - 8<sup>1</sup>/<sub>2</sub>" (16cm - 22cm)  
 6<sup>7</sup>/<sub>8</sub>" - 9<sup>1</sup>/<sub>4</sub>" (18cm - 24cm)

#### Bella Strong Gauntlet

Size

Palm



S 5<sup>7</sup>/<sub>8</sub>" - 7" (15cm - 18cm)  
 M 7" - 8<sup>1</sup>/<sub>4</sub>" (18cm - 21cm)  
 L 8<sup>1</sup>/<sub>4</sub>" - 9<sup>1</sup>/<sub>2</sub>" (21cm - 24cm)  
 XL 9<sup>1</sup>/<sub>2</sub>" - 10<sup>1</sup>/<sub>2</sub>" (24cm - 27cm)

Wrist



5" - 6<sup>3</sup>/<sub>4</sub>" (13cm - 17cm)  
 5<sup>1</sup>/<sub>2</sub>" - 7<sup>1</sup>/<sub>2</sub>" (14cm - 19cm)  
 6<sup>1</sup>/<sub>4</sub>" - 8<sup>1</sup>/<sub>4</sub>" (16cm - 21cm)  
 7" - 9" (18cm - 23cm)

Material	Class II 20-30 mmHg	Colors Available
Elvarex Seamless Glove	<input type="checkbox"/>	Beige <input type="checkbox"/>

**The Non-Custom Elvarex Seamless 20-30 mmHg Glove with finger stubs is a less expensive alternative for patients paying out of pocket. If your patient has hand and finger swelling, you should consider getting a custom-made glove. \*Please note that custom-made gloves are not seamless.**

### Measuring and ordering information:

Size	Beige	Compression**	Palm	Wrist
JOBST Elvarex Seamless Glove				
1	7955600	20 - 30mmHg	6 <sup>1</sup> / <sub>8</sub> "-6 <sup>7</sup> / <sub>8</sub> " (16cm-17.5cm)	5"-5 <sup>3</sup> / <sub>4</sub> " (13cm-14.5cm)
2	7955601	20 - 30mmHg	6 <sup>7</sup> / <sub>8</sub> "-7 <sup>1</sup> / <sub>2</sub> " (17.5cm-19cm)	5 <sup>3</sup> / <sub>4</sub> "-6 <sup>1</sup> / <sub>4</sub> " (14.5cm-16cm)
3	7955602	20 - 30mmHg	7 <sup>1</sup> / <sub>2</sub> "-8" (19cm-20.5cm)	6 <sup>1</sup> / <sub>4</sub> "-6 <sup>7</sup> / <sub>8</sub> " (16cm-17.5cm)
4	7955603	20 - 30mmHg	8"-8 <sup>5</sup> / <sub>8</sub> " (20.5cm-22cm)	6 <sup>7</sup> / <sub>8</sub> "-7 <sup>1</sup> / <sub>2</sub> " (17.5cm-19cm)
5	7955604	20 - 30mmHg	8 <sup>5</sup> / <sub>8</sub> "-9 <sup>1</sup> / <sub>4</sub> " (22cm-23.5cm)	7 <sup>1</sup> / <sub>2</sub> "-8" (19cm-20.5cm)
6	7955605	20 - 30mmHg	9 <sup>1</sup> / <sub>4</sub> "-9 <sup>3</sup> / <sub>4</sub> " (23.5cm-25cm)	8"-8 <sup>5</sup> / <sub>8</sub> " (20.5cm-22cm)



Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1360 N. Sandburg Terrace · Suite #103 · Chicago, IL 60610 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com