



Account Number: 1014233

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

JUZO NON-CUSTOM CIRCULAR KNIT ARM SLEEVE MEASUREMENT FORM AND SIZING CHART

Measure in centimeters

Left Right

Length Measurement from axilla to wrist (c-g) = _____

LONG LENGTH IS > 43 CMS. REGULAR LENGTH IS < 43 CMS.

The length measurement is taken from **point 'c'** along the inside of the arm to the inseam, **point 'g'**, of the arm; measure with patient standing or sitting with arm down by the side.

Measurement Points

Circumference **point 'g'** is taken around the upper arm at the axilla. For ease in determining this point, place a book up into the arm pit area. Point 'g' on the arm will be even with the top of the book.


Circumference **point 'e'** is found in the middle of the elbow.

Circumference **point 'c'** is the point of greatest compression and therefore the most important point. This point is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point. **Be sure not to measure too tight.**

Circumference 'g' = _____ **c-g with full silicone** * *Soft 2000 only available with full silicone.*
Circumference 'e' = _____ **c-g with no silicone**
Circumference 'c' = _____ **c-h over the shoulder**

Material	Class I 20-30 mmHg	Class II 30-40 mmHg	Colors Available
Dynamic (Varin Soft-In)	3511 <input type="checkbox"/>	3512 <input type="checkbox"/>	Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/>
Soft 2000	2001 <input type="checkbox"/>	2002 <input type="checkbox"/>	Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/> Cinnamon <input type="checkbox"/> _____ <input type="checkbox"/> * Seasonal colors available

* To view seasonal colors, visit www.dreamsleeve.com



C-G/C-H Regular - C-G length <17" / 43 cm
C-G/C-H Long - C-G length >17" / 43 cm

	I	II	III	IV	V	VI
g max	26.5-33 cm 10 1/2"-13"	28.5-35.0 cm 11 1/4"-13 3/4"	30.5-37.5 cm 12"-14 3/4"	32.5-39.5 cm 12 3/4"-15 1/2"	34.0-42.0 cm 13 1/2"-16 1/2"	36.0-44.0 cm 14 1/2"-17 1/2"
g	21.5-26.5 cm 8 1/2"-10 1/2"	23.5-28.5 cm 9 1/4"-11 1/4"	25.0-31.0 cm 9 3/4"-12 1/4"	27.0-33.0 cm 10 1/2"-13"	29.0-35.0 cm 11 1/2"-13 3/4"	30.5-37.5 cm 12"-14 3/4"
e max	22.5-27.5 cm 8 3/4"-10 3/4"	24.5-29.5 cm 9 3/4"-11 1/2"	26-32 cm 10 1/4"-12 1/2"	27.5-33.5cm 10 3/4"-13 1/4"	29.0-35.0 cm 11 1/2"-13 3/4"	30.0-37.0 cm 11 3/4"-14 1/2"
e	20.5-25.5 cm 8"-10"	22-27cm 8 3/4"-10 3/4"	23.5-28.5 cm 9 1/4"-11 1/4"	25.0-30.0 cm 9 3/4"-11 3/4"	26.0-32.0cm 10 1/4"-12 1/2"	27.5-33.5 cm 10 3/4"-13 1/4"
c	14.0-15.5 cm 5 1/2"-6"	15.5-17.0 cm 6"-6 3/4"	17.0-18.5 cm 6 3/4"-7 1/4"	18.5-20 cm 7 1/4"-7 3/4"	20.0-21.5 cm 7 3/4"-8 1/2"	21.5-23 cm 8 1/2"-9"