



Account Number: 1014233

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

**Measure in Centimeters**

**JUZO NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASUREMENT FORM AND SIZING CHART**

Left  Right

Thigh Length \_\_\_\_\_

Calf Length \_\_\_\_\_

A-g REGULAR LENGTH IS 72-83 CM  
 A-g SHORT LENGTH 63-72 CM  
 A-g PETITE LENGTH 55-63 CM

A-d REGULAR LENGTH IS 40-46 CM  
 A-d SHORT LENGTH IS 34-40 CM  
 A-d PETITE LENGTH IS <33 CM

**For thigh length**, the length measurement is taken from the floor (at the inside of heel) along the inside of the leg to the inseam/groin, point 'g'; this length should end at the gluteal fold.

**For calf length**, the length measurement is taken from the floor (at the inside of heel) along the inside of the leg to approximately 3 1/2 centimeters below the crease of the knee, point 'd'

Silicone borders are not usually recommended for calf length stockings, however, they shaped legs.

**For calf and thigh length stockings**, measure with the patient standing.

**Measurement Points**

Circumference **point 'g'** is found at the inseam/groin

Circumference **point 'c'** is the largest circumference of the calf

Circumference **point 'b'** is the point of greatest compression and therefore the most important point. This point is the smallest circumference of the ankle found closest to and above the ankle bone

**Circumference 'g'** = \_\_\_\_\_

**Circumference 'c'** = \_\_\_\_\_

**Circumference 'b'** = \_\_\_\_\_

Open Toe  Silicone

Closed Toe  No Silicone


**Measurement Point for the Body Part (KM) or Hip Attachment**

\*The following measurement will be necessary if you are ordering a "custom-made" thigh length stocking to be sewn to a body part.

\*It is common for a Class 1 Body Part (KM) to be sewn to a Class 2 or 3 leg stocking.

**Hip circumference** = \_\_\_\_\_

KM Class 1 (18-21mmHg)  KM Class 2 (23-32mmHg)



Material <small>*Dynamic is new name for Varin</small>	Class 1 20-30 mmHg	Class 2 30-40 mmHg	Class 3 40-50 mmHg	Colors Available
Dynamic (calf length) AD	3511 <input type="checkbox"/>	3512 <input type="checkbox"/>	3513 <input type="checkbox"/>	Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/>
Dynamic (thigh length) AG	3511 <input type="checkbox"/>	3512 <input type="checkbox"/>	3513 <input type="checkbox"/>	Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/>
Dynamic Full Knit (thigh length) AG		3512L <input type="checkbox"/>	3513L <input type="checkbox"/>	Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/>
Dynamic Full Knit (panty hose)			3513ATLY <input type="checkbox"/>	
Soft 2000 (calf length) AD	2001 <input type="checkbox"/>	2002 <input type="checkbox"/>		Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/> _____ <input type="checkbox"/> * Seasonal colors available
Soft 2000 (thigh length) AG	2001 <input type="checkbox"/>	2002 <input type="checkbox"/>		Beige <input type="checkbox"/> Chestnut <input type="checkbox"/> Black <input type="checkbox"/> _____ <input type="checkbox"/> * Seasonal colors available

\* To view seasonal colors, visit [www.juzousa.com](http://www.juzousa.com), then click on Product > Stocking > Fashion > Soft 2000