



Account Number: 1014233

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

Measure in centimeters

Left Right

JUZO NON-CUSTOM HAND GLOVE/GAUNTLET MEASUREMENT FORM AND SIZING CHART

Measurement Points:

Circumference **point 'a'** is found across the metacarpus.
 Circumference **point 'c'** is found distally on the arm where the wrist meets the hand.
 One should find the smallest circumference at this point.

Circumference **'a'** = _____

Circumference **'c'** = _____

Style	
Gauntlet w/ thumb stub only	<input type="checkbox"/>
Glove w/ finger stubs	<input type="checkbox"/>

Material	Class 1 18-21 mmHg	Class 2 23-32 mmHg	Color Available
Expert *formerly Helastic	3021 <input type="checkbox"/>	3022 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/>


Non-custom made gloves with finger stubs cost the same for us to provide as the **custom-made gloves with finger stubs**.
 Finger lengths end at the PIP joints for non-custom gloves. If your patient has hand and finger swelling, you should consider a custom-made glove whereby the finger lengths terminate at the DIP joints.

Material	Class I 20-30 mmHg	Class II 30-40 mmHg	Colors Available
1100 Series Gauntlet	1101 <input type="checkbox"/>	1102 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/>
Dreamsleeve Soft 2000 Gauntlet	2001 <input type="checkbox"/>	2002 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/> Cinnamon <input type="checkbox"/> <small>_____ <input type="checkbox"/> * Seasonal colors available</small>

* 1100 Series and Soft 2000 are only available as gauntlet with thumb stub.

* To view seasonal colors, visit www.dreamsleeve.com

Size Chart for Juzo® 3021/3022 Gauntlet



	1		2		3		4		5		6	
	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm
①	6.5	5	6.5	5	6.5	5	7.5	6	7.5	6	7.5	6
②	6.2	5	6.7	5.5	7.1	5.9	7.6	6.4	8	6.8	8.5	7.3
③	6	5	6.3	5.3	6.7	5.7	7.1	6	7.5	6.4	7.6	6.7
④	5.7	4	6	4.2	6.3	4.5	6.7	4.8	7.1	5.1	7.4	5.4
⑤	5.5	4	5.9	4.2	6.3	4.5	6.6	4.8	7	5.1	7.3	5.4
a	17.5 cm 6 7/8"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"		23.5 cm 9 1/4"		25 cm 9 3/4"	
c	14.5 cm 5 3/4"		16 cm 6 1/4"		17.5 cm 7"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"	

Size Chart for Juzo® 1101 / 1102, 2001 / 2002 Gauntlet

	1-S	2-M	3-L	4-XL
a	17.5-19.5 cm 6 7/8"-7 5/8"	19.5-21.0 cm 7 5/8"-8 1/4"	21.0-23.0 cm 8 1/4"-9 1/8"	23.0-25.5 cm 9 1/8"-10 1/8"
c	14.5-16.5 cm 5 5/8"-6 1/2"	16.5-18.0 cm 6 1/2"-7"	18.0-20.0 cm 7"-7 7/8"	20.0-22 cm 7 7/8"-8 5/8"