



Account Number: 120167

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____



MEDI NON-CUSTOM CIRCULAR KNIT & FLAT KNIT HAND GLOVE MEASUREMENT FORM AND SIZING CHART

Measure in centimeters

Left Right

Measurement Points (Two measurements are required)

Circumference **point 'a'** is found across the metacarpus (hand).

Circumference **point 'c'** is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point (wrist).

Circumference **'a'** = _____

Circumference **'c'** = _____

*Note that the measurements above will provide us the information to order a **hand glove with finger stubs**. The finger lengths will extend to the DIP joints.

Item	Class I	Class II	Color(s) Available
Medi 95	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	Beige <input type="checkbox"/> Mocha <input type="checkbox"/>
Medi 95	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	Beige <input type="checkbox"/>
Medi 95	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	Caramel <input type="checkbox"/>

* Medi 95 & Medi have half compression in fingers; no seams.

* Medi Espirit has full compression in fingers; with seams.

SIZE CHART

Medi 95 & Medi

	II	III	IV	V	VI	VII
A	16-19 cm 6¼" - 7½"	18-21 cm. 7" - 8¼"	20-23 cm. 8" - 9"	22-25 cm. 8¾" - 9¾"	24- 27cm. 9½" - 10¾"	26-29 cm. 10¼" - 11½"
C	13-15 cm 5" - 6"	15-17 cm. 6" - 6¾"	17-19 cm. 6¾" - 7½"	19-21 cm. 7½" - 8¼"	21-23 cm. 8¼" - 9"	23-25 cm. 9" - 9¾"

Flat Knit - Medi Espirit

	I	II	III	IV	V
A	17.5-18.5 cm 6¾" - 7¼"	18.5 - 20 cm 7¼" - 7¾"	20-22 cm 7¾" - 8¾"	22-24 cm 8¾" - 9½"	24- 27cm 9½" - 10¾"
C	15-17 cm 5¾" - 6¾"	17-19 cm 6¾" - 7½"	19-21 cm 7½" - 8¼"	21-23 cm. 8¼" - 9¼"	23-25 cm 9¼" - 10¼"

