



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

MEDI NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASUREMENT FORM AND SIZING CHART



Measure in Centimeters Left Right

Thigh Length _____ Calf Length _____

A-g STANDARD LENGTH IS > 71 CM
 A-g PETITE LENGTH < 71 CM

A-d STANDARD LENGTH IS > 38 CM
 A-d PETITE LENGTH IS < 38 CM

For thigh length, the length measurement is taken from the floor (inside of heel) along the inside of the leg to the inseam/groin, point 'g'; this length should end at the gluteal fold.


For calf length, the length measurement is taken from the floor (inside of heel) along the inside of the leg to approximately 3 1/2 centimeters below the crease of the knee, point 'd'

Measurement Points

Circumference point 'g' is found at the inseam/groin

Circumference point 'c' is the largest circumference of the calf

Circumference point 'b' is the point of greatest compression and therefore the most important point. This point is the smallest circumference of the ankle found closest to and above the ankle bone



Circumference 'g' = _____

Measurement Point for Waist Attachment

The following measurement will be necessary if you are ordering a thigh length stocking to be sewn to a waist attachment

Circumference 'c' = _____

Hip Circumference 't' = _____

Waist attachment option not available in Medi Forte stockings.

Circumference 'b' = _____

Item	Class I	Class II	Class III	Options**
Medi Plus	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	40-50 mmHg <input type="checkbox"/>	Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe* <input type="checkbox"/>
Medi Forte	N/A	30-40 mmHg <input type="checkbox"/>	40-50 mmHg <input type="checkbox"/>	Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/> N/A

*Closed toe option not available in Medi Plus, Class III or Medi Forte stockings.

**Medi Plus and Medi Forte stockings only come in beige.

	I	II	III	IV	V	VI	VII
g panty & thigh	38-57 cm 15"-22½"	41-62 cm 16"-24¼"	44-67 cm. 17¼"-26¼"	47-72 cm. 18½"- 28¼"	50-77 cm. 19¾"-30¼"	53-82 cm. 21"-32¼"	56-87 cm. 22¼"- 34¼"
g thigh w/ silicone only	42-49 cm 16½"-19¼"	45-53 cm 17¾"- 20¾"	48-57 cm. 19"-22¼"	52-61 cm. 20½"- 24"	56-65 cm. 22"- 25¾"	60-70 cm. 23½"- 27½"	64-74 cm. 25"- 29"
c calf & calf w/ silicone	26-34 cm 10¼"-13¼"	28-37 cm 11"-14½"	30-40 cm. 11¾"-15¾"	32-43 cm. 12½"- 17"	34-46 cm. 13¼"-18"	36-49 cm. 14¼"- 19¼"	38-52 cm. 15"- 20½"
c *calf extra wide & calf extra wide w/ silicone	34-38 cm 13¼"-15"	37-42 cm 14½"- 16½"	40-46 cm. 15¾"-18"	43-50 cm. 17"- 19½"	46-54 cm. 18"-21"	49-58 cm. 19¼"- 22¾"	52-62 cm. 20½"- 24½"
b ankle	16.5-18.5 cm 6½"-7¼"	19-21 cm 7½"-8¼"	21.5-23.5 cm. 8½"-9¼"	24-26 cm. 9½"- 10¼"	26.5-28.5 cm. 10½"-11¼"	29-32 cm. 11½"- 12½"	32.5-35 cm. 12¾"- 13¾"

*Medi Plus extra wide calf-length stockings are **only available in classes I and II.**

Medi Forte extra wide calf-length stockings are **only available in classes II and III.