



Account Number: 120167

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Fitter Title: _____ (example PT/OT/PTA)

Date: _____

MEDI NON-CUSTOM CIRCULAR KNIT ARM SLEEVE MEASUREMENT FORM AND SIZING CHART

Measure in centimeters

Left Right

Length Measurement from axilla to wrist (c-g) = _____



LONG LENGTH IS > 44 CM.

STANDARD LENGTH IS < 44 CM.

The length measurement is taken from **point 'c'** along the inside of the arm to the inseam, **point 'g'**, of the arm; measure with patient sitting and arm extended straight on table.

Measurement Points

Circumference **point 'g'** is taken around the upper arm at the axilla. For ease in determining this point, place a folded piece of paper up into the arm pit area. Point 'g' on the arm will be even with the top of the paper. **Measure with moderate tension.**

Circumference **point 'e'** is found in the middle of the elbow. **Measure without tension.**

Circumference **point 'c'** is the point of greatest compression and therefore the most important point. This point is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point. **Measure without tension.**


Circumference 'g' = _____ **c-g with silicone border**

Circumference 'e' = _____ **c-g with no silicone**

Circumference 'c' = _____

Item	Class I	Class II	Colors Available	Note
Medi 95	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	Beige <input type="checkbox"/> Mocha <input type="checkbox"/>	If no silicone is selected, sleeve only comes in beige.

SIZE CHART



	II	III	IV	V	VI	VII
G	22-29 cm	25-32 cm	29-36 cm	33-40 cm	37-44 cm	41-48 cm
E	20-26 cm	23-29 cm	26-32 cm	29-35 cm	32-38 cm	35-41 cm
C	13-15 cm	15-17 cm	17-19 cm	19-21 cm	21-23 cm	23-25 cm