



Account Number: 120167

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Fitter Title: _____ (example PT/OT/PTA)

Date: _____

MEDI NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASUREMENT FORM AND SIZING CHART



Measure in Centimeters Left Right

Thigh Length _____

Calf Length _____

A-g STANDARD LENGTH IS > 71 CM
A-g PETITE LENGTH < 71 CM

A-d STANDARD LENGTH IS > 38 CM
A-d PETITE LENGTH IS < 38 CM

For thigh length, the length measurement is taken from the floor (inside of heel) along the inside of the leg to the inseam/groin, point 'g'; this length should end at the gluteal fold.


For calf length, the length measurement is taken from the floor (inside of heel) along the inside of the leg to approximately 3 1/2 centimeters below the crease of the knee, point 'd'

Measurement Points

Circumference **point 'g'** is found at the inseam/groin

Circumference **point 'c'** is the largest circumference of the calf

Circumference **point 'b'** is the point of greatest compression and therefore the most important point. This point is the smallest circumference of the ankle found closest to and above the ankle bone



Circumference 'g' = _____

Measurement Point for Waist Attachment

The following measurement will be necessary if you are ordering a thigh length stocking to be sewn to a waist attachment

Circumference 'c' = _____

Hip Circumference 't' = _____

Waist attachment option not available in Medi Forte stockings.

Circumference 'b' = _____

Item	Class I	Class II	Class III	Options**
Medi Plus	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	40-50 mmHg <input type="checkbox"/>	Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe* <input type="checkbox"/>
Medi Forte	N/A	30-40 mmHg <input type="checkbox"/>	40-50 mmHg <input type="checkbox"/>	Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/> N/A

*Closed toe option not available in Medi Plus, Class III or Medi Forte stockings.

**Medi Plus and Medi Forte stockings only come in beige.

	I	II	III	IV	V	VI	VII
g panty & thigh	38-57 cm 15"-22½"	41-62 cm 16"-24¼"	44-67 cm. 17¼"-26¼"	47-72 cm. 18½"- 28¼"	50-77 cm. 19¾"-30¼"	53-82 cm. 21"-32¼"	56-87 cm. 22¼"- 34¼"
g thigh w/ silicone only	42-49 cm 16½"-19¼"	45-53 cm 17¾"- 20¾"	48-57 cm. 19"-22¼"	52-61 cm. 20½"- 24"	56-65 cm. 22"- 25¾"	60-70 cm. 23½"- 27½"	64-74 cm. 25"- 29"
c calf & calf w/ silicone	26-34 cm 10¼"-13¼"	28-37 cm 11"-14½"	30-40 cm. 11¾"-15¾"	32-43 cm. 12½"- 17"	34-46 cm. 13¼"-18"	36-49 cm. 14¼"- 19¼"	38-52 cm. 15"- 20½"
c *calf extra wide & calf extra wide w/ silicone	34-38 cm 13¼"-15"	37-42 cm 14½"- 16½"	40-46 cm. 15¾"-18"	43-50 cm. 17"- 19½"	46-54 cm. 18"-21"	49-58 cm. 19¼"- 22¾"	52-62 cm. 20½"- 24½"
b ankle	16.5-18.5 cm 6½"-7¼"	19-21 cm 7½"-8¼"	21.5-23.5 cm. 8½"-9¼"	24-26 cm. 9½"- 10¼"	26.5-28.5 cm. 10½"-11¼"	29-32 cm. 11½"- 12½"	32.5-35 cm. 12¾"- 13¾"

*Medi Plus extra wide calf-length stockings are **only available in classes I and II.**

Medi Forte extra wide calf-length stockings are **only available in classes II and III.

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1360 N. Sandburg Terrace · Suite #103 · Chicago, IL 60610 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com