



SIGVARIS Prescription

PATIENT INFORMATION

Patient	Date	Number of Pairs	
Diagnosis	Duration of Treatment	Number of Refills	<input type="checkbox"/> Substitutions Allowed


Armsleeve

w/ Gauntlet
 w/out Gauntlet

Compression Level	Accessories
<input type="checkbox"/> 15–20mmHg	<input type="checkbox"/> Latex-Free Gloves
<input type="checkbox"/> 20–30mmHg	<input type="checkbox"/> Rubber Gloves
<input type="checkbox"/> 30–40mmHg	<input type="checkbox"/> Cornu-Thenard Extensor (Donning Device)
	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)

Physician Name	License Number
Physician Signature X	Phone
	Fitter

UPPER EXTREMITY



ARMSLEEVE WITHOUT GAUNTLET

- 1 Wrist Circumference at narrowest part of wrist (above wrist bone) **(c)** _____
- 2 Axilla (underarm) Circumference just underneath the armpit. **(f)** _____
- 3 Arm length from wrist to axilla (measure outside of arm only) **(C-G)** _____

ARMSLEEVE WITH GAUNTLET

(continue from steps 1–3)

- 4 Palm Circumference at the widest part of hand (thumb excluded) **(a)** _____

Send to:

Luna Medical, Inc.

1816 W. Belmont Ave., Ste. 1
Chicago, IL 60657

P: 800.380.4339 F: 888.696.0299

lunamedical.com

ADVANCE ARMSLEEVE SIZING CHART

PALM (a)	SMALL 6.5"–8" (17–21cm)		MEDIUM 7"–8.5" (18–22cm)				LARGE 7.5"–9" (19–23cm)				EXTRA-LARGE 8"–9.5" (20–24cm)			
WRIST (c)	SMALL 5.5"–6.5" (14–16cm)		MEDIUM 6.5"–7" (16–18cm)				LARGE 7"–8" (18–20cm)				EXTRA-LARGE 8"–9" (20–23cm)			
UPPER MID-ARM (f)	REGULAR 10"–13.5" (26–34cm)		REGULAR 11.5–14.5" (29–37cm)		PLUS 14.5"–17" (37–43cm)		REGULAR 12.5–15.5" (32–40cm)		PLUS 15.5"–18" (40–46cm)		REGULAR 13.5–16.5" (34–42cm)		PLUS 16.5"–19" (42–48cm)	
LENGTH (C–G)	REGULAR ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)	REG. ≤17" (≤43cm)	LONG >17" (>43cm)
SIZE	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

For indications and contraindications, please see prescribing information at sigvarisusa.com

SIGVARIS is a registered trademark of SIGVARIS AG, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2015 SIGVARIS, Inc.



SIGVARIS Prescription

PATIENT INFORMATION


Patient	Date	Number of Pairs	
Diagnosis	Duration of Treatment	Number of Refills	<input type="checkbox"/> Substitutions Allowed

Style				
<input type="checkbox"/> Closed Toe		<input type="checkbox"/> Open Toe		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/>
Calf	Thigh	Pantyhose	Thigh w/ waist	Maternity/ Plus Sizes

Compression Level	Accessories
<input type="checkbox"/> 15–20mmHg	<input type="checkbox"/> Latex-Free Gloves
<input type="checkbox"/> 18–25mmHg	<input type="checkbox"/> Rubber Gloves
<input type="checkbox"/> 20–30mmHg	<input type="checkbox"/> Cotton Underliner
<input type="checkbox"/> 30–40mmHg	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)
<input type="checkbox"/> 40–50mmHg	<input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)
<input type="checkbox"/> 50–60mmHg	<input type="checkbox"/> Cornu-Thenard Extensor (Donning Device)

Physician Name	License Number
Physician Signature X	Phone
	Fitter

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone (**b**) _____
- 2 Calf Circumference at fullest part of the calf (**c**) _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) (**A-D**) _____

PANTYHOSE AND THIGH-HIGH STYLES

(continue from steps 1–3)

- 4 Thigh Circumference at its widest circumference (**g**) _____
- 5 Thigh Length from floor up to the gluteal fold (**A-G**) _____
- 6 Hip Circumference at widest part of the hip (**f**) _____

Send to:

Luna Medical, Inc.

1816 W. Belmont Ave., Ste. 1
Chicago, IL 60657

P: 800.380.4339 F: 888.696.0299
lunamedical.com

WOMEN'S SHOE SIZES:

PETITE LEG 4–8, SHORT LEG 5.5–9.5, LONG LEG 8–12

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14

PLUS SIZES HIP CIRCUMFERENCE GUIDELINES:

SMALL: <56" (142cm), MEDIUM: <62" (157.5cm),
LARGE: <70" (178cm)

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)			MEDIUM 8.5"–10" (21.5–25.5cm)			LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)			MEDIUM 12.5"–17.5" (31.5–45cm)			LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	PETITE* <14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	PETITE* <14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18–26" (45.5–66cm)			MEDIUM 19–28.5" (48–73cm)			LARGE 20.5–33" (52–84cm)		*SP & MP SIZES ARE AVAILABLE IN CLOSED TOE ONLY, WITHOUT GRIP-TOP	
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5–44.5" (83–112.5cm)			MEDIUM 34.5–50" (87–127cm)			LARGE 38.5–59.5" (98–150.5cm)			
LEG LENGTH: WOMEN & OPEN TOE	PETITE* ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5 (72cm)	PETITE* ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5 (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)		
LEG LENGTH: MEN	SHORT ≤30" (76cm)		LONG >30 (76cm)	SHORT ≤30" (76cm)		LONG >30 (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)		
SIZE	SP*	SS	SL	MP*	MS	ML	LS	LL	XS	XL

For indications and contraindications, please see prescribing information at sigvarisusa.com

SIGVARIS is a registered trademark of SIGVARIS AG, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2015 SIGVARIS, Inc.