



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

SIGVARIS

* Please check off options *

Style				
<input type="checkbox"/> Closed Toe	<input type="checkbox"/> Open Toe			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/>
Calf	Thigh	Pantyhose	Thigh w/ waist	Maternity/ Plus Sizes

Men's Shoe Size _____ Women's Shoe Size _____

WOMEN'S SHOE SIZES:


PETITE LEG 4-8, SHORT LEG 5.5-9.5, LONG LEG 8-12

MEN'S SHOE SIZES: SHORT LEG 7-11, LONG LEG 9.5-14

PLUS SIZES HIP CIRCUMFERENCE GUIDELINES:

SMALL: <56" (142CM), MEDIUM: <62" (157.5CM), LARGE: <70" (178CM)

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone (b) _____
- 2 Calf Circumference at fullest part of the calf (c) _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) (A-D) _____

PANTYHOSE AND THIGH-HIGH STYLES

(continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference (g) _____
- 5 Thigh Length from floor up to the gluteal fold (A-G) _____
- 6 Hip Circumference at widest part of the hip (f) _____

Material	15-20 mHg	20-30 mHg	30-40 mHg	40-50 mHg	Colors
Sheer Fashion For Women	120 <input type="checkbox"/>	N/A	N/A	N/A	<input type="checkbox"/> Charcoal <input type="checkbox"/> Natural <input type="checkbox"/> Cafe <input type="checkbox"/> Taupe <input type="checkbox"/> Suntan <input type="checkbox"/> Black
Business Casual For Men	189 <input type="checkbox"/>	N/A	N/A	N/A	<input type="checkbox"/> Navy <input type="checkbox"/> Charcoal <input type="checkbox"/> Brown <input type="checkbox"/> Black
Casual Cotton	186 Mens <input type="checkbox"/> 146 Womens <input type="checkbox"/>	N/A	N/A	N/A	<input type="checkbox"/> White <input type="checkbox"/> Brown <input type="checkbox"/> Navy <input type="checkbox"/> Khaki <input type="checkbox"/> Black
Cotton	N/A	232 Mens <input type="checkbox"/> 232 Womens <input type="checkbox"/>	233 Mens <input type="checkbox"/> 233 Womens <input type="checkbox"/>	N/A	<input type="checkbox"/> Navy <input type="checkbox"/> Black <input type="checkbox"/> Chocolate <input type="checkbox"/> White <input type="checkbox"/> Crispa
Cushioned Cotton	182 Mens <input type="checkbox"/> 142 Womens <input type="checkbox"/>	362 Mens <input type="checkbox"/> 362 Womens <input type="checkbox"/>	N/A	N/A	<input type="checkbox"/> Black <input type="checkbox"/> White
Sea Island Cotton	151 Mens <input type="checkbox"/> 141 Womens <input type="checkbox"/>	222 Mens <input type="checkbox"/> 222 Womens <input type="checkbox"/>	N/A	N/A	<input type="checkbox"/> Black <input type="checkbox"/> Brown <input type="checkbox"/> Navy
All Season Motrino Wool	192 Mens <input type="checkbox"/> 152 Womens <input type="checkbox"/>	242 Mens <input type="checkbox"/> 242 Womens <input type="checkbox"/>	N/A	N/A	<input type="checkbox"/> Navy <input type="checkbox"/> Brown <input type="checkbox"/> Black
Allure	711 <input type="checkbox"/>	712 <input type="checkbox"/>	N/A	N/A	<input type="checkbox"/> Black <input type="checkbox"/> Graphite
Eversheer	781 <input type="checkbox"/>	782 <input type="checkbox"/>	783 <input type="checkbox"/>	N/A	<input type="checkbox"/> Dark Navy <input type="checkbox"/> Suntan <input type="checkbox"/> Black <input type="checkbox"/> Natural <input type="checkbox"/> Nightshade
Soft Opaque	841 <input type="checkbox"/>	842 <input type="checkbox"/>	843 <input type="checkbox"/>	N/A	<input type="checkbox"/> Midnight Blue <input type="checkbox"/> Pecan <input type="checkbox"/> Black <input type="checkbox"/> Nude <input type="checkbox"/> Graphite <input type="checkbox"/> Espresso <input type="checkbox"/> Natural <input type="checkbox"/> Crispa <input type="checkbox"/> Black * 20-30 mHg only <input type="checkbox"/> Suntan <input type="checkbox"/> Black
Select Comfort	N/A	862 <input type="checkbox"/>	863 <input type="checkbox"/>	N/A	<input type="checkbox"/> Dark Navy <input type="checkbox"/> Pink <input type="checkbox"/> Onyx <input type="checkbox"/> Stripe <input type="checkbox"/> Stripe
Microfiber Shades	N/A	832 Mens <input type="checkbox"/> 832 Womens <input type="checkbox"/>	N/A	N/A	<input type="checkbox"/> Tan Khaki <input type="checkbox"/> Black <input type="checkbox"/> Steel Grey * Men's Only * Calf Only
Midtown Microfiber For Men	821 <input type="checkbox"/>	822 <input type="checkbox"/>	823 <input type="checkbox"/>	N/A	
Access	N/A	922 <input type="checkbox"/> 972 <input type="checkbox"/>	923 <input type="checkbox"/> 973 <input type="checkbox"/>	N/A	<input type="checkbox"/> Beige <input type="checkbox"/> Brown
Natural Rubber	N/A	N/A	503 <input type="checkbox"/>	504 <input type="checkbox"/>	Beige Only / Opentoe Only

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1057 W. Grand Ave · Suite 1 · Chicago, IL 60642 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com