




Account Number: 4057807

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____



ReadyWrapTM
 MANUFACTURED BY SOLARIS

LOWER EXTREMITY SIZING FORM



NEW

NEW

NEW

ReadyWrap Thigh		Length		Girths
Item Code	Item Desc.	Thigh	Above Knee	Groin
RW-LE-EG-S25	Average, Small	up to 25 cm	38 - 48 cm	48 - 63 cm
RW-LE-EG-M25	Average, Medium	up to 25 cm	45 - 58 cm	58 - 73 cm
RW-LE-EG-L25	Average, Large	up to 25 cm	52 - 68 cm	68 - 83 cm
RW-LE-EG-XL25	Average, X-Large	up to 25 cm	59 - 78 cm	78 - 93 cm
RW-LE-EG-XXL25	Average, XX-Large	up to 25 cm	66 - 88 cm	88 - 103 cm
RW-LE-EG-S30	Tall, Small	up to 30 cm	38 - 48 cm	48 - 63 cm
RW-LE-EG-M30	Tall, Medium	up to 30 cm	45 - 58 cm	58 - 73 cm
RW-LE-EG-L30	Tall, Large	up to 30 cm	52 - 68 cm	68 - 83 cm
RW-LE-EG-XL30	Tall, X-Large	up to 30 cm	59 - 78 cm	78 - 93 cm
RW-LE-EG-XXL30	Tall, XX-Large	up to 30 cm	66 - 88 cm	88 - 103 cm

ReadyWrap Knee		Girths	
Item Code	Item Desc.	Below Knee	Above Knee
RW-LE-DE-S	Small	28 - 38 cm	38 - 48 cm
RW-LE-DE-M	Medium	33 - 43 cm	45 - 58 cm
RW-LE-DE-L	Large	38 - 48 cm	52 - 68 cm
RW-LE-DE-XL	X-Large	43 - 58 cm	59 - 78 cm
RW-LE-DE-XXL	XX-Large	53 - 68 cm	66 - 88 cm

ReadyWrap Calf		Length		Girths
Item Code	Item Desc.	Calf	Ankle	Calf
RW-LE-BD-S30	Average, Small	up to 30 cm	18 - 23 cm	28 - 38 cm
RW-LE-BD-M30	Average, Medium	up to 30 cm	20 - 28 cm	33 - 43 cm
RW-LE-BD-L30	Average, Large	up to 30 cm	25 - 33 cm	38 - 48 cm
RW-LE-BD-XL30	Average, X-Large	up to 30 cm	30 - 38 cm	43 - 58 cm
RW-LE-BD-XXL30	Average, XX-Large	up to 30 cm	35 - 43 cm	53 - 68 cm
RW-LE-BD-S35	Tall, Small	up to 35 cm	18 - 23 cm	28 - 38 cm
RW-LE-BD-M35	Tall, Medium	up to 35 cm	20 - 28 cm	33 - 43 cm
RW-LE-BD-L35	Tall, Large	up to 35 cm	25 - 33 cm	38 - 48 cm
RW-LE-BD-XL35	Tall, X-Large	up to 35 cm	30 - 38 cm	43 - 58 cm
RW-LE-BD-XXL35	Tall, XX-Large	up to 35 cm	35 - 43 cm	53 - 68 cm

ReadyWrap Foot		Length		Girths
Item Code	Item Desc.	Foot	Y	Ankle
RW-LE-AB-L/R-S	Left / Right, Small	up to 21 cm	28 - 33 cm	18 - 23 cm
RW-LE-AB-L/R-M	Left / Right, Medium	up to 23 cm	30 - 38 cm	20 - 28 cm
RW-LE-AB-L/R-L	Left / Right, Large	up to 26 cm	35 - 43 cm	25 - 33 cm
RW-LE-AB-L/R-XL	Left / Right, X-Large	up to 29 cm	40 - 48 cm	30 - 38 cm

ReadyWrap Toe		Girth
Item Code	Item Desc.	MTP
RW-LE-AA-L/R-S	Left / Right, Small	18.0 - 20.4 cm
RW-LE-AA-L/R-M	Left / Right, Medium	20.5 - 23.5 cm
RW-LE-AA-L/R-L	Left / Right, Large	23.6 - 26.0 cm
RW-LE-AA-L/R-XL	Left / Right, X-Large	26.1 - 28.5 cm
RW-LE-AA-L/R-XXL	Left / Right, XX-Large	28.6 - 31.0 cm