




## UPPER EXTREMITY SIZING FORM

All sizes in centimeters



READYWRAP ARM			SIZING		
Item Code	Item Desc.	Length	Wrist	Elbow	Axilla
RW-UE-CG-S43	Short Small	43	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M43	Short Medium	43	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L43	Short Large	43	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL43	Short X-Large	43	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0
RW-UE-CG-S46	Average Small	46	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M46	Average Medium	46	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L46	Average Large	46	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL46	Average X-Large	46	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0
RW-UE-CG-S49	Tall Small	49	13.5 - 17.0	19.0 - 26.0	22.0 - 29.0
RW-UE-CG-M49	Tall Medium	49	15.5 - 20.0	24.0 - 33.0	26.0 - 36.0
RW-UE-CG-L49	Tall Large	49	17.5 - 23.0	29.0 - 40.0	31.0 - 43.0
RW-UE-CG-XL49	Tall X-Large	49	17.5 - 23.0	30.0 - 43.0	36.0 - 50.0

READYWRAP GAUNTLET			SIZING	
Item Code	Item Desc.	Length	Palm	Wrist
RW-UE-BC-L/R-S	Left/Right Small	13	18.0 - 20.0	13.5 - 17.0
RW-UE-BC-L/R-M	Left/Right Med.	14	20.0 - 22.0	15.5 - 20.0
RW-UE-BC-L/R-L	Left/Right Large	15	22.0 - 24.0	17.5 - 23.0
RW-UE-BC-L/R-XL	Left/Right X-Large	16	24.0 - 26.0	19.5 - 26.0

Patient Last Name: \_\_\_\_\_

Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_

Fitter First Name: \_\_\_\_\_

Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)

Date: \_\_\_\_\_

Item Code	Item Description	Black Qty.	Beige Qty.	Item Code	Item Description	Black Qty.	Beige Qty.
<b>ReadyWrap Arm</b>				<b>ReadyWrap Gauntlet</b>			
RW-UE-CG-S43	Short Small			RW-UE-BC-L-S	Left Small		
RW-UE-CG-M43	Short Med.			RW-UE-BC-R-S	Right Small		
RW-UE-CG-L43	Short Large			RW-UE-BC-L-M	Left Med.		
RW-UE-CG-XL43	Short XL			RW-UE-BC-R-M	Right Med.		
RW-UE-CG-S46	Avg. Small			RW-UE-BC-L-L	Left Large		
RW-UE-CG-M46	Avg. Med.			RW-UE-BC-R-L	Right Large		
RW-UE-CG-L46	Avg. Large			RW-UE-BC-L-XL	Left XL		
RW-UE-CG-XL46	Avg. XL			RW-UE-BC-R-XL	Right XL		
RW-UE-CG-S49	Tall Small						
RW-UE-CG-M49	Tall Med.						
RW-UE-CG-L49	Tall Large						
RW-UE-CG-XL49	Tall XL						