



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____




CIRCAID JUXTAFIT FOOT OPTIONS (AFW) MEASUREMENT FORM

FOOT OPTIONS

*afw = ankle foot wrap

	Small	Medium/Standard	Large	X-Large
compression anklets				
customizable interlocking afw*				
juxtalite afw*				
juxtafit premium afw*				
juxtafit premium interlocking afw*				
pac band™				
single band afw*				



circaid power added compression band (pac band)

circumferences in cm.			
	standard	large	xlarge
cB	19 - 27	25 - 50	<56
cA1	19 - 36	22 - 41	<44


circaid compression anklet

circumferences in cm.	
	large
cB	25 - 42
cA1	22 - 41


*circaid compression anklet only available in standard and large.

circaid® juxtafit™ premium ankle foot wrap (afw)

circumferences and lengths in cm.			
size	S	M	L
ℓ A - Y	> 9	> 9	> 9
cB	< 36	< 36	< 41
cH	19 - 24	24 - 29	29 - 34
ℓ Z	> 16	> 16	> 16



circaid® juxtalite™ ankle foot wrap (afw)




circaid® juxtalite™ ankle foot wrap (afw)

circumferences in cm.		
	S	M
cH	19 - 26	24 - 31

NOTE: The maximum foot arch circumference for the juxtalite™ ankle foot wrap is 31cm. If your patient exceeds this measurement, the single band ankle foot wrap™ is available for order (see opposing page).

circaid® single band ankle foot wrap (afw)

dimensions in cm.	
width	7.6
length	203



circaid juxtafit premium interlocking ankle foot wrap (afw)

circumferences and lengths in cm.			
size	S	M	L
ℓ A - Y	> 9	> 9	> 9
cB	< 36	< 36	< 41
cA1	19 - 24	24 - 29	29 - 34
ℓ Z	13 - 25	13 - 25	13 - 25