



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_



# CIRCAID NON-CUSTOM JUXTAFIT ESSENTIALS LEG GARMENT MEASUREMENT FORM

## juxtafit essentials lower leg

	S	M	M-FC	L	L-FC	XL	XL-FC	XXL
Short (28 cm)								
Long (36 cm)								

## juxtafit essentials upper leg with knee

	XXS		XS		S		M		L		XL	
	L	R	L	R	L	R	L	R	L	R	L	R
X-Short (35 cm)												
Short (45 cm)												
Long (55 cm)												

## capri

A	B	C	D	E	F	G	H	I	J

## FOOT OPTIONS

\*afw = ankle foot wrap

	Small	Medium/Standard	Large	X-Large
compression anklets				
customizable interlocking afw*				
juxtalite afw*				
juxtafit premium afw*				
juxtafit premium interlocking afw*				
pac band™				
single band afw*				

### circaid juxtafit premium & essentials lower leg

circumference in centimeters								
size	S	M	M-X full calf	L	L-X full calf	XL	XL-X full calf	XXL
cC	30-38	36-43	43-50	41-48	48-56	41-51	51-61	53-64
cB1	20-28	25-33	33-40	30-38	38-46	33-41	41-51	43-53
cB	15-23	20-28	20-28	25-33	25-33	30-38	30-38	33-41


  

actual product length in centimeters		
size	short	long
cB-D	28	36

leg length		
size	short	long
cB-D	<36	>36

All products contain the circaid Built-In-Pressure (BPS) guide card. Product lengths may differ from body measurements. Size small juxtafit receive standard compression ankle. All other sizes receive large compression ankle. Size small, medium, medium-full calf, and large juxtafit receive small cover up. All other sizes receive large cover up. Size small, medium, and medium-full calf juxtafit receive standard pac band. All other sizes receive large pac band.



### circaid juxtafit essentials upper leg with knee

circumference in centimeters						
size	XXS	XS	S	M	L	XL
cG	40-52	46-62	54-70	60-80	69-90	77-101
cE1	32-42	37-48	42-55	50-65	58-76	66-87
cD	23-30	27-35	31-41	35-46	40-52	44-58

actual product length in centimeters			
size	x-short	short	long
cD-G	35	45	55

sizing measurement			
size	x-short	short	long
cE-G	<25	<35	>35

Product lengths may differ from body measurements.

### circaid capri

size	circumference in inches		circumference in centimeters	
	waist (T)	thigh (G)	waist (T)	thigh (G)
A	25-31	18-24	63-79	46-62
B	26-33	19-26	67-85	49-65
C	28-36	20-27	72-91	52-70
D	30-38	22-29	76-96	55-73
E	32-40	24-32	81-102	60-81
F	33-43	25-33	85-108	63-84
G	35-44	26-35	90-113	67-90
H	37-47	27-37	94-119	70-93
I	39-49	28-38	99-125	72-96
J	41-51	30-40	103-130	77-102

When taking measurements, measure at the largest circumference point. When selecting a size, choose the smallest size match for the most snug fit. Choose the largest size match for a more relaxed, easy-to-wear fit.