



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____



CIRCAID NON-CUSTOM JUXTALITE LEG GARMENT MEASUREMENT FORM

LOWER EXTREMITY


juxtalite® lower leg

	S	M	M-FC	L	L-FC	XL	XL-FC	XXL
Short (28 cm)								
Long (33 cm)								

FOOT OPTIONS

	Small	Medium/Standard	Large	X-Large
juxtalite afw*				

*afw = ankle foot wrap




circaid juxtalite lower leg

circumference in centimeters								
size	S	M	M-X	L	L-X	XL	XL-X	XXL
cC	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64
cB	19-22	22-27	22-27	25-32	25-32	30-37	30-37	35-42
actual product length in centimeters								
size	short				long			
floor-D	28				33			
sizing measurement	short			long				
floor-D1	<44			>44				

Product lengths may differ from body measurements.

circaid® juxtalite™ ankle foot wrap (afw)



circaid® juxtalite™ ankle foot wrap (afw)

circumferences in cm.		
	S	M
cH	19-26	24-31



NOTE: The maximum foot arch circumference for the juxtalite™ ankle foot wrap is 33cm. If your patient exceeds this measurement, the single band ankle foot wrap™ is available for order (see opposing page).