



Account Number: 120167  
 Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

**CIRCAID** CUSTOM ARMSLEEVE SIZING CHART  
*Making Compression Easier.*

The Juxta-Fit™ Armsleeve is constructed from the wrist to the axilla measurement. A lateral rise is incorporated at the top of the garment in addition to a 10° bend is at the elbow. Measure patient's arm length with the arm bent at 35°. Also ensure that the arm is at its smallest and maintain arm size at time of measuring to assure a proper fit of the garment(s). **NOTE: 1in = 2.54cm**

ARM AND HAND SELECTION	
<input type="checkbox"/> RIGHT	<input type="checkbox"/> LEFT

COVERUP COLOR
<input type="checkbox"/> BLACK
<input type="checkbox"/> BEIGE
<input type="checkbox"/> PINK

JUXTA-FIT™ GAUNTLET SIZE OPTIONS					
PALM CIRCUMFERENCE	14-17cm	17-20cm	20-23cm	23-26cm	26-29cm
SIZE	<input type="checkbox"/> XSMALL	<input type="checkbox"/> SMALL	<input type="checkbox"/> MEDIUM	<input type="checkbox"/> LARGE	<input type="checkbox"/> XL

ARM LENGTH	
c-e	c-g
wrist to elbow crease	wrist to axilla (or where the garment should end)

HAND OPTIONS
<input type="checkbox"/> GAUNTLET
<input type="checkbox"/> GLOVE
<input type="checkbox"/> DORSUM STRAP
<input type="checkbox"/> NONE

JUXTA-FIT™ GLOVE MEASUREMENTS											
LENGTH	1	2	3	4	5	6	7	8	9	10	
	1x	2x	3x	4x	5x	a	b	c			
CIRCUMFERENCE	1w	2w	3w	4w	5w						

ARM CIRCUMFERENCE MEASUREMENTS (Palm; Wrist; 5cm increments from the wrist to the axilla)	MEASURING POINTS
	Palm
	Wrist
	5cm
	10cm
	15cm
	20cm
	25cm
	30cm
	35cm
	40cm
	45cm
	50cm
	55cm

