



Patient Last Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_  
 Date: \_\_\_\_\_

Account Number: 120167  
 Patient First Name: \_\_\_\_\_  
 Fitter First Name: \_\_\_\_\_  
 (example PT/OT/PTA)



## Non-Custom Lower Leg Garment Measurement Form & Sizing Chart

### SIZING CHART

Patient should be standing when taking measurements. Limbs should be at their smallest; limb size should be maintained from time of measurements to ensure a proper fit of garment. When on borderline between sizes, select a smaller size if significant edema is present and reduction is anticipated; if edema is minimal, select larger size. After measuring the smallest dimension of the ankle and the largest dimension of the calf, if any of the measurements do not fall within one of these size ranges, call (800) 247-2243. Length measurement is for Juxta-Fit™ only, actual legging lengths are 28cm (Short) and 36cm (Long); Ready-Fit™ and Classic-Flex™ have one adjustable length. **Note 1in = 2.54cm**

MEASUREMENTS		
	LEFT	RIGHT
<b>B</b>		
<b>B<sub>1</sub></b>		
<b>C</b>		
<b>H</b>		
<b>B-D</b>		

JUXTA-FIT™ PREMIUM & ESSENTIALS			
SIZE	CIRCUMFERENCE MEASUREMENT		
	B	B <sub>1</sub>	C
<b>SMALL</b>	15-23cm	20-28cm	30-38cm
<b>MEDIUM</b>	20-28cm	25-33cm	36-43cm
<b>MEDIUM FULL CALF</b>	20-28cm	33-40cm	43-50cm
<b>LARGE</b>	25-33cm	30-38cm	41-48cm
<b>LARGE FULL CALF</b>	25-33cm	38-46cm	48-56cm
<b>XL</b>	30-38cm	33-41cm	41-51cm
<b>XL FULL CALF</b>	30-38cm	41-51cm	51-61cm
<b>XXL*</b>	33-41cm	43-53cm	53-64cm

AFW MEASUREMENT	
SIZE	H
<b>SMALL</b>	19-24cm
<b>MEDIUM</b>	24-29cm
<b>LARGE</b>	29-34cm


LENGTH MEASUREMENT	
SIZE	B-D
<b>SHORT</b>	under 36cm
<b>LONG</b>	over 36cm

LEGGING OPTIONS	SIZE:	QTY:
<b>28CM JUXTA-FIT™ PREMIUM † STANDARD LOWER LEGGING (SHORT)</b>		
<b>36CM JUXTA-FIT™ PREMIUM † STANDARD LOWER LEGGING (LONG)</b>		

ANCILLARIES	SIZE:	QTY:
COMFORT COVERUP™ FOR LOWER LEG (BLACK) (Small: 34cm - 50cm)		
COMFORT COVERUP™ FOR LOWER LEG (BLACK) (Large: 51cm - 80cm)		
COMFORT COVERUP™ FOR LOWER LEG (BEIGE) (Small: 34cm - 50cm)		
COMFORT COVERUP™ FOR LOWER LEG (BEIGE) (Large: 51cm - 80cm)		
COMFORT CAPRI™		
(For a sizing chart visit <a href="http://www.circaid.com">www.circaid.com</a> )		
INSTRUCTIONAL DVD		

\* If measurements fall within XXL size range, please consider ordering a custom Juxta-Fit™ for optimal fit.  
 † EZ Single-Band Ankle-Foot Wrap™ included by default  
 ‡ Juxta-Lite™ Ankle-Foot Wrap is available only in sizes Small and Medium  
 † Comfort Compression Ankle™ included by default

FOOT OPTIONS	SIZE:	QTY:
Comfort Accessory Kit (3 Pairs of Standard Comfort Compression Ankle™, 3 Pairs of Large Comfort Compression Ankle™, 3 Comfort Leg Liners™ & 3 Comfort PAC Bands™)		
Comfort Compression Ankle™ (Standard max circumference H = 36cm B = 27cm) (Large max circumference H = 41cm B = 50cm)		
Comfort PAC Band™ (Standard max circumference H = 36cm B = 27cm) (Large max circumference H = 41cm B = 50cm)		
Comfort EZ Single-Band Ankle-Foot Wrap™ (3")	N/A	
Comfort EZ Single-Band Ankle-Foot Wrap™ (4")	N/A	
Standard Juxta-Fit™ Ankle-Foot Wrap (max H circumference 34cm)		
Standard Juxta-Fit™ Interlocking Ankle-Foot Wrap (max H circ. 34cm)		
Custom Juxta-Fit™ Interlocking Ankle-Foot Wrap		
Juxta-Lite™ Ankle-Foot Wrap (max H circumference 31cm) ‡		



PRICES SUBJECT TO CHANGE WITHOUT NOTICE.