




**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

## circaid® reduction kit sizing information for leg



circaid® reduction kit - upper leg			
Circumference in cm	regular	wide	
<b>cG</b>	<80	<125	
<b>cE2</b>	<65	<110	
Length in cm	short	standard	long
<b>l KC - Groin</b>	<37	37-41	>41

### Each upper leg kit includes:

- 1 reduction upper leg component
- 1 reduction shelf strap
- 6 Velcro® stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

circaid® reduction kit - lower leg			
Circumference in cm	regular	wide	
<b>cC</b>	<65	<80	
<b>cB</b>	<55	<65	
Length in cm	short	standard	long
<b>l mall - KC</b>	<32	32-36	>36

### Each lower leg kit includes:

- 1 reduction lower leg component
- 1 reduction kit undersleeve leg
- 1 reduction shelf strap
- 6 Velcro® stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

circaid® reduction kit - knee	
Circumference in cm	regular
<b>cE2</b>	<110
<b>cE</b>	<102
<b>cC</b>	<80
one length available:	standard 30cm

### Each knee kit includes:

- 1 reduction knee component
- 2 Velcro® knee spines
- 1 paper measuring tape
- 1 direction for use
- 1 reduction kit shelf strap

circaid® reduction kit - foot options				
Circumference in cm	compression anklet		pac band*	single band AFW
	standard	large	x-large	
<b>cH</b>	19-36	22-41	<44	>44
<b>cB</b>	19-27	25-50	<56	>56

**Accessories list:** • reduction kit undersleeve leg • reduction kit shelf strap • undersleeve arm standard and wide