




# FarrowWrap® 4000



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

### 3 Easy Measuring Steps:

1. Length (B-D) - Ankle to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) - Calf (widest)
3. Circumference (B) - Ankle (smallest part)



Measurements		
	Left	Right
cC		
cB		
/B-D		

### JOBST® FarrowWrap® 4000

Ref No.	Size	Length	Color	Quantity
7666100	XS	Regular	Tan	
7666101	S	Regular	Tan	
7666102	M	Regular	Tan	
7666103	L	Regular	Tan	
7666108	XS	Regular	Black	
7666109	S	Regular	Black	
7666110	M	Regular	Black	
7666111	L	Regular	Black	
7666104	XS	Tall	Tan	
7666105	S	Tall	Tan	
7666106	M	Tall	Tan	
7666107	L	Tall	Tan	
7666112	XS	Tall	Black	
7666113	S	Tall	Black	
7666114	M	Tall	Black	
7666115	L	Tall	Black	

### JOBST® FarrowWrap® 4000 Size Chart


	XS	S	M	L
<b>C</b> Mid-Calf	25-38 cm	32-48 cm	37-56 cm	43-67 cm
<b>B</b> Ankle	17-23 cm	19-28 cm	25-37 cm	29-43 cm
<b>B-D</b> Regular	30-33 cm	30-33 cm	30-33 cm	30-33 cm
<b>B-D</b> Tall	34-37 cm	34-37 cm	34-37 cm	34-37 cm



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

### 3 Easy Measuring Steps:

1. Length (A-D) - Floor to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) - Calf (widest)
3. Circumference (B) - Ankle (smallest part)



Measurements		
	Left	Right
cC		
cB		
∅A-D		

### JOBST® FarrowWrap® Basic (30-40 mmHg\*)

Ref No.	Size	Length	Quantity
7666000	XS	Regular	
7666001	S	Regular	
7666002	M	Regular	
7666003	L	Regular	
7666004	XL	Regular	
7666005	XS	Tall	
7666006	S	Tall	
7666007	M	Tall	
7666008	L	Tall	
7666009	XL	Tall	

### JOBST® FarrowWrap® Basic Size Chart

	XS	S	M	L	XL
<b>C</b> Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
<b>B</b> Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
<b>A-D</b> Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
<b>A-D</b> Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm

\* The mean compression for an average size ankle  
 Please refer to the product label and / or package insert for full instructions on the safe use of these products.