




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ (example: PT/OT/PTA)
Date: _____

FarrowWrap® Liner/ Farrowhybrid® Liner Order Form



JOBST FARROWWRAP TG® SOFT KNEE - HIGH LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666900	<40 cm (C)	TG SOFT AD (LEG) SIZE SMALL	
7666902	40-70 cm (C)	TG SOFT AD (LEG) SIZE MEDIUM	
7666904	70-120 cm (C)	TG SOFT AD (LEG) SIZE LARGE	

JOBST FARROWWRAP TG® SOFT THIGH - HIGH LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7667000	<40 cm (G)	TG SOFT AG (THIGH) SIZE SMALL	
7667002	40-70 cm (G)	TG SOFT AG (THIGH) SIZE MEDIUM	
7667004	70-120 cm (G)	TG SOFT AG (THIGH) SIZE LARGE	

JOBST® FARROWFOAM® LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666500	X SMALL	SIDE CLOSURE AD LINERS	
7666501	SMALL	SIDE CLOSURE AD LINERS	
7666502	MEDIUM	SIDE CLOSURE AD LINERS	
7666503	LARGE	SIDE CLOSURE AD LINERS	
7666504	X LARGE	SIDE CLOSURE AD LINERS	
7666505	XX LARGE	SIDE CLOSURE AD LINERS	

JOBST FARROWFOAM® LINERS SIZING CHART

SIZE	ANKLE (CB)	CALF (CC)	FIBULAR HEAD (CD)	LENGTH (A-D)
XS	18-23 cm	32-38 cm	37-43 cm	38 cm
S	21-26 cm	37-43 cm	42-48 cm	39 cm
M	23-28 cm	42-48 cm	45-51 cm	41 cm
L	26-31 cm	47-53 cm	50-56 cm	42 cm
XL	28-33 cm	52-58 cm	55-61 cm	43 cm
XXL	31-36 cm	57-63 cm	60-66 cm	45 cm

FARROWHYBRID® SIZING CHART

A' _____	Standard:	20-27cm
	Wide:	28-35cm

JOBST® FARROW LINER

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666800	15-58 cm (C)	SILVER LINER AD (LEG) SIZE SMALL	
7666802	52-70 cm (C)	SILVER LINER AD (LEG) SIZE LARGE	



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



FarrowWrap® Liner/ FarrowHybrid® Liner Order Form

JOBST® FARROWHYBRID® FOOT COMPRESSION FARROWLINER™

BNR CODE	DESCRIPTION	SIZE	LENGTH	COLOR	SIZING	QTY
7666400	FarrowHybrid® ADI Foot Compression 20-30mmHg	Small	Standard	Taupe	mens shoe size 3-6 womens shoe size 4-7	
7666401	FarrowHybrid® ADI Foot Compression 20-30mmHg	Medium	Standard	Taupe	mens shoe size 7-9 womens shoe size 8-10	
7666402	FarrowHybrid® ADI Foot Compression 20-30mmHg	Large	Standard	Taupe	mens shoe size 10-12 womens shoe size 11-13	
7666403	FarrowHybrid® ADI Foot Compression 20-30mmHg	Extra Large	Standard	Taupe	mens shoe size 13-15 womens shoe size 14-16	
7666412	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Small	Wide	Taupe	mens shoe size 3-6 womens shoe size 4-7	
7666413	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Medium	Wide	Taupe	mens shoe size 7-9 wide womens shoe size 8-10 wide	
7666414	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Large	Wide	Taupe	mens shoe size 10-12 wide womens shoe size 11-13 wide	
7666415	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Extra Large	Wide	Taupe	mens shoe size 13-15 wide womens shoe size 14-16 wide	
7666416	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Small	Wide	Black	mens shoe size 3-6 wide womens shoe size 4-7 wide	
7666417	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Medium	Wide	Black	mens shoe size 7-9 wide womens shoe size 8-10 wide	
7666418	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Large	Wide	Black	mens shoe size 10-12 wide womens shoe size 11-13 wide	
7666419	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Extra Large	Wide	Black	mens shoe size 13-15 wide womens shoe size 14-16 wide	
7666404	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Small	Standard	Sand	mens shoe size 3-6 womens shoe size 4-7	
7666405	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Medium	Standard	Sand	mens shoe size 7-9 womens shoe size 8-10	
7666406	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Large	Standard	Sand	mens shoe size 10-12 womens shoe size 11-13	
7666407	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Extra Large	Standard	Sand	mens shoe size 13-15 womens shoe size 14-16	
7666408	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Small	Wide	Sand	mens shoe size 3-6 wide womens shoe size 4-7 wide	
7666409	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Medium	Wide	Sand	mens shoe size 7-9 wide womens shoe size 8-10 wide	
7666410	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Large	Wide	Sand	mens shoe size 10-12 wide womens shoe size 11-13 wide	
7666411	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Extra Large	Wide	Sand	mens shoe size 13-15 wide womens shoe size 14-16 wide	

PLEASE NOTE: All liners available in Taupe unless otherwise noted in the above size description.

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies
 1057 W. Grand Ave · Suite 1 · Chicago, IL 60642 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com