




# FarrowWrap® Ready-to-Wear



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_



CIRCUMFERENCE (C)	LEFT	RIGHT
cG (groin)		
cE (above patella)		
cD (2 finger widths below crease)		
cC (widest calf)		
cB (smallest point above ankle)		
cA <sup>1</sup> (widest foot)		





LENGTH	LEFT	RIGHT
/E-G (medially)		
/A-D (posteriorly, follow contour)		
/X (from base of toes to heel)		

Base of small toe to heel.

**Please note:**

- FarrowWrap® Basic comes with FarrowHybrid liner
- All other legpieces come with silver liners
- Footpiece comes with silver liner
- Thighpiece comes with TG soft liner
  - Classic & Strong Thighpiece come with Strong Kneepiece
  - Lite Thighpiece comes with Lite Kneepiece
- Extra liners can be ordered for additional charge on the Liner order form

READY-TO-WEAR					
Thighpiece	XS	S	M	L	XL
<b>G</b> Groin	60-70 cm	65-75 cm	70-80 cm	75-90 cm	85-100 cm
<b>E</b> Above Knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm
<b>E-G</b> Short	16-20 cm	16-20 cm	16-20 cm	16-20 cm	16-20 cm
<b>E-G</b> Regular	21-26 cm	21-26 cm	21-26 cm	21-26 cm	21-26 cm
<b>E-G</b> Tall	27-31 cm	27-31 cm	27-31 cm	27-31 cm	27-31 cm
Legpiece	XS	S	M	L	XL
<b>C</b> Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
<b>B</b> Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
<b>A-D</b> Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
<b>A-D</b> Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm
Footpiece	XS	S	M	L	XL
<b>A<sup>1</sup></b> Midfoot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
<b>X</b> Regular	16-17 cm	17-18 cm	19-20 cm	20-21 cm	22-23 cm
<b>X</b> Tall	18-19 cm	19-20 cm	21-22 cm	22-23 cm	24-25 cm
Kneepiece	XS	S	M	L	XL
<b>E</b> Above Knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm

<b>JOBST® FARROWWRAP® LITE (20-30 mmHg)</b>				
Footpiece	BSN Code	Size	Length	Qty.
	7665205	XS	Reg	
	7665206	S	Reg	
	7665207	M	Reg	
	7665208	L	Reg	
	7665209	XL	Reg	
	7665200	XS	Long	
	7665201	S	Long	
	7665202	M	Long	
	7665203	L	Long	
	7665204	XL	Long	
Legpiece	BSN Code	Size	Length	Qty.
	7665400	XS	Reg	
	7665401	S	Reg	
	7665402	M	Reg	
	7665403	L	Reg	
	7665404	XL	Reg	
	7665405	XS	Tall	
	7665406	S	Tall	
	7665407	M	Tall	
	7665408	L	Tall	
	7665409	XL	Tall	
Thighpiece	BSN Code	Size	Length	Qty.
	7665505	XS	Short	
	7665506	S	Short	
	7665507	M	Short	
	7665508	L	Short	
	7665509	XL	Short	
	7665500	XS	Reg	
	7665501	S	Reg	
	7665502	M	Reg	
	7665503	L	Reg	
	7665504	XL	Reg	
	7665510	XS	Tall	
	7665511	S	Tall	
	7665512	M	Tall	
	7665513	L	Tall	
	7665514	XL	Tall	
Kneepiece	BSN Code	Size	Length	Qty.
	7336700	XS	N/A	
	7336701	S	N/A	
	7336702	M	N/A	
	7336703	L	N/A	
	7336704	XL	N/A	

<b>JOBST® FARROWWRAP® STRONG (30-40 mmHg)</b>			
BSN Code	Size	Length	Qty.
7664604	XS	Reg	
7664605	S	Reg	
7664606	M	Reg	
7664607	L	Reg	
7664608	XL	Reg	
7664609	XS	Long	
7664600	S	Long	
7664601	M	Long	
7664602	L	Long	
7664603	XL	Long	
BSN Code	Size	Length	Qty.
7664800	XS	Reg	
7664801	S	Reg	
7664802	M	Reg	
7664803	L	Reg	
7664804	XL	Reg	
7664805	XS	Tall	
7664806	S	Tall	
7664807	M	Tall	
7664808	L	Tall	
7664809	XL	Tall	
BSN Code	Size	Length	Qty.
7664905	XS	Short	
7664906	S	Short	
7664907	M	Short	
7664908	L	Short	
7664909	XL	Short	
7664900	XS	Reg	
7664901	S	Reg	
7664902	M	Reg	
7664903	L	Reg	
7664904	XL	Reg	
7664910	XS	Tall	
7664911	S	Tall	
7664912	M	Tall	
7664913	L	Tall	
7664914	XL	Tall	
BSN Code	Size	Length	Qty.
7336800	XS	N/A	
7336801	S	N/A	
7336802	M	N/A	
7336803	L	N/A	
7336804	XL	N/A	

<b>JOBST® FARROWWRAP® CLASSIC (30-40 mmHg)</b>			
BSN Code	Size	Length	Qty.
7665605	XS	Reg	
7665606	S	Reg	
7665607	M	Reg	
7665608	L	Reg	
7665609	XL	Reg	
7665600	XS	Long	
7665601	S	Long	
7665602	M	Long	
7665603	L	Long	
7665604	XL	Long	
BSN Code	Size	Length	Qty.
7665700	XS	Reg	
7665701	S	Reg	
7665702	M	Reg	
7665703	L	Reg	
7665704	XL	Reg	
7665705	XS	Tall	
7665706	S	Tall	
7665707	M	Tall	
7665708	L	Tall	
7665709	XL	Tall	
BSN Code	Size	Length	Qty.
7665806	XS	Short	
7665807	S	Short	
7665808	M	Short	
7665809	L	Short	
7665810	XL	Short	
7665800	XS	Reg	
7665801	S	Reg	
7665802	M	Reg	
7665803	L	Reg	
7665804	XL	Reg	
7665812	XS	Tall	
7665813	S	Tall	
7665814	M	Tall	
7665815	L	Tall	
7665816	XL	Tall	

<b>JOBST® FARROWWRAP® BASIC (30-40 mmHg)</b>			
BSN Code	Size	Length	Qty.
7665905	XS	Reg	
7665906	S	Reg	
7665907	M	Reg	
7665908	L	Reg	
7665909	XL	Reg	
7665900	XS	Tall	
7665901	S	Tall	
7665902	M	Tall	
7665903	L	Tall	
7665904	XL	Tall	
BSN Code	Size	Length	Qty.
7666000	XS	Reg	
7666001	S	Reg	
7666002	M	Reg	
7666003	L	Reg	
7666004	XL	Reg	
7666005	XS	Tall	
7666006	S	Tall	
7666007	M	Tall	
7666008	L	Tall	
7666009	XL	Tall	

**Please note:**  
**This side of form**  
**must be submitted**  
**with front side.**