




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

FarrowWrap® Trim-To-Fit Knee High Order Form



FarrowWrap® TTF Knee High

FarrowWrap® TTF Legpiece		Medium	XLarge	FarrowWrap® TTF Footpiece		Medium
(L4) & (L5) Circumferences		20 – 60 cm	45 – 90 cm	Forefoot (X) circumference		20 – 30 cm
(L3) Circumference		20 – 55 cm	45 – 90 cm	(X) Heel to forefoot length		14 – 21 cm
(L2) Circumference		20 – 45 cm	35 – 80 cm	LITE TTF Footpiece		FWLT-T-FT3 Qty
(B) Ankle Circumference		20 – 40 cm	35 – 70 cm	STRONG TTF Footpiece		FWST-T-FT3 Qty
2B-D) Posterior Leg Length		27 – 35 cm	27 – 35 cm	Please note that you may substitute a FarrowHybrid® Foot Compression for a footpiece. See FarrowHybrid® order form for additional details.		
LITE TTF Legpiece	FWLT-T-LT3	Qty				
STRONG TTF Legpiece	FWST-T-LT	Qty	Qty			

FarrowWrap® TTF garments fit a very wide range of sizes and can be trimmed on site for a perfect fit, and retrimmed if the patient's edema further reduces.

When measuring the circumferences from B to L5, the distance between each circumference will vary for the different fabrics. For **LITE** legpieces, measure the circumference at **point B** and then measure the circumference every **6.5 cm** until measurement **point D** is reached. For **STRONG** legpieces, measure the circumference at **point B** and then measure the circumference every **5 cm** until measurement **point D** is reached/

FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. When a TTF footpiece and legpiece are purchased together, a free pair of liners is included. Silver liners are supplied by default, but you may request Terry Cloth instead.

THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY