



FarrowWrap® Upper Extremity

Ready-To-Wear • Custom • Trim-to-Fit



Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Fitter Title: _____ (example: PT/OT/PTA)

Date: _____

Ready-To-Wear Arm

Fabric

LITE

Color

Beige (00)

Black (01)

Length

Short (S)

Regular (R)

Long (L)

Size

Small (2)

Medium (3)

Large (4)

Side

Left (L)

Right (R)

Custom

Style

Custom

Fabric

Classic

LITE

STRONG

Color

Tan

Side

Left (L)

Right (R)

Trim-To-Fit

Size


Medium (3)

Extra Large (5)

Side

Left (L)

Right (R)



Circumference (c)	Left	Right
cG (axilla)*		
cF ¹ (widest bicep)		
cF (distal bicep)		
cE (elbow crease)*		
cD ¹ (widest forearm)		
cD (distal forearm)		
cC (wrist)*		

Lengths	Left	Right
IC-G (wrist to axilla)*		
IC-E (wrist to elbow)		

Comments

*Circumference measurements needed for Ready to Wear garments. Please note that all measurements are needed for Custom garments.

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1057 W. Grand Ave · Suite 1 · Chicago, IL 60642 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com



FarrowWrap® Upper Extremity

Ready-To-Wear • Custom • Trim-to-Fit

READY-TO-WEAR ARM			
Size	Wrist (cC)	Elbow (cE)	Axilla (cG)
S	14-18 cm	20-27 cm	22-31 cm
M	16-21 cm	25-34 cm	29-39 cm
L	19-25 cm	30-40 cm	32-45 cm
Length (C-G)	Short 40-43 cm		
	Regular 44-47 cm		
	Long 48-51 cm		

TRIM-TO-FIT ARM							
Size	Wrist (cC)	Forearm (cD)	Elbow (cE)	Bicep (cF)	Axilla (cG)	Length (C-E)	Length (E-G)
M	15-31 cm	20-36 cm	20-38 cm	25-42 cm	25-44 cm	20-25 cm	15-23 cm
XL	15-35 cm	20-40 cm	25-42 cm	25-45 cm	30-48 cm	20-28 cm	

READY-TO-WEAR ARM - TAN		
BNR	Description	Quantity
7665035	LITE Armpiece Long Small Left Tan	
7665000	LITE Armpiece Long Small Right Tan	
7665001	LITE Armpiece Long Medium Left Tan	
7665002	LITE Armpiece Long Medium Right Tan	
7665003	LITE Armpiece Long Large Left Tan	
7665004	LITE Armpiece Long Large Right Tan	
7665005	LITE Armpiece Regular Tan Small Left	
7665006	LITE Armpiece Regular Tan Small Right	
7665007	LITE Armpiece Regular Tan Medium Left	
7665008	LITE Armpiece Regular Tan Medium Right	
7665009	LITE Armpiece Regular Tan Large Left	
7665010	LITE Armpiece Regular Tan Large Right	
7665011	LITE Armpiece Short Tan Small Left	
7665012	LITE Armpiece Short Tan Small Right	
7665013	LITE Armpiece Short Tan Medium Left	
7665014	LITE Armpiece Short Tan Medium Right	
7665015	LITE Armpiece Short Tan Large Left	
7665016	LITE Armpiece Short Tan Large Right	

READY-TO-WEAR ARM - BLACK		
BNR	Description	Quantity
7665023	LITE Armpiece Regular Black Small Left	
7665024	LITE Armpiece Regular Black Small Right	
7665025	LITE Armpiece Regular Black Medium Left	
7665026	LITE Armpiece Regular Black Medium Right	
7665027	LITE Armpiece Regular Black Large Left	
7665028	LITE Armpiece Regular Black Large Right	
7665029	LITE Armpiece Short Black Small Left	
7665030	LITE Armpiece Short Black Small Right	
7665031	LITE Armpiece Short Black Medium Left	
7665032	LITE Armpiece Short Black Medium Right	
7665033	LITE Armpiece Short Black Large Left	
7665034	LITE Armpiece Short Black Large Right	
7665017	LITE Armpiece Long Black Small Left	
7665018	LITE Armpiece Long Black Small Right	
7665019	LITE Armpiece Long Black Medium Left	
7665020	LITE Armpiece Long Black Medium Right	
7665021	LITE Armpiece Long Black Large Left	
7665022	LITE Armpiece Long Black Large Right	

CUSTOM CLASSIC		
BNR	Description	Quantity
7599500	Classic Arm	

CUSTOM STRONG		
BNR	Description	Quantity
7598700	STRONG Arm	

CUSTOM LITE		
BNR	Description	Quantity
7599100	LITE Custom Arm (Left or Right)	

FARROWWRAP TRIM-TO-FIT		
BNR	Description	Quantity
7612000	LITE Armpiece (Left) Medium with 1 pair of silver arm liners	
7612002	LITE Armpiece (Right) Medium with 1 pair of silver arm liners	
7612001	LITE Armpiece (Left) Extra Large with 1 pair of silver arm liners	
7612003	LITE Armpiece (Right) Extra Large with 1 pair of silver arm liners	