




# FarrowWrap® Custom



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_



| CIRCUMFERENCE (C)                 | LEFT | RIGHT |
|-----------------------------------|------|-------|
| cG (groin)                        |      |       |
| cF (mid thigh)                    |      |       |
| cE (above patella)                |      |       |
| cD (2 finger widths below crease) |      |       |
| cC (widest calf)                  |      |       |
| cB² (midway between cB¹ and cC)   |      |       |
| cB¹ (base of calf)                |      |       |
| cB (smallest point above ankle)   |      |       |
| cA¹ (widest foot)                 |      |       |

| LENGTH                             | LEFT | RIGHT |
|------------------------------------|------|-------|
| lE-G (medially)                    |      |       |
| lA-D (posteriorly, follow contour) |      |       |
| lX (from base of toes to heel)     |      |       |

Base of small toe to heel.

| JOBST® FARROWWRAP® LINERS |   |          |
|---------------------------|---|----------|
| BSN Code                  | Description   | Quantity |
| 7666700                   | Farrow Silver AG Liners (Arm Hand Circumference < 56cm)   |          |
| 7666600                   | Farrow Silver CG Liners (Arm Circumference < 56cm)        |          |
| 7666800                   | Farrow Silver AD Liners Small (Leg Circumference < 56cm)  |          |
| 7666802                   | Farrow Silver AD Liners Large (Leg Circumference 53-70cm) |          |
| 7666900                   | Farrow TG® Soft AD Liners (S < 40cm)                      |          |
| 7666902                   | Farrow TG® Soft AD Liners (M 40-70cm)                     |          |
| 7666904                   | Farrow TG® Soft AD Liners (L > 70cm)                      |          |
| 7667000                   | Farrow TG® Soft AG Liners (S < 40cm)                      |          |
| 7667002                   | Farrow TG® Soft AG Liners (M 40-70cm)                     |          |
| 7667004                   | Farrow TG® Soft AG Liners (L > 70cm)                      |          |
| 7667006                   | Farrow TG® Soft Arm Liners (S < 40cm)                     |          |
| 7667008                   | Farrow TG® Soft Arm Liners (M 40-70cm)                    |          |
| 7667010                   | Farrow TG® Soft Arm Liners (L > 70cm)                     |          |

| BRAND / COMPRESSION                          | QUANTITY  | ADDITIONAL LINERS (ordered in pairs)  |
|--|---|---|
| <input type="checkbox"/> LITE (20-30mmHg)    | Thigh: Left _____ Right _____   | NOTE: there is a charge for these liners.<br><br>Liner model:<br><input type="checkbox"/> TG Soft <input type="checkbox"/> Silver<br><br>Size:<br><input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large |
| <input type="checkbox"/> CLASSIC (30-40mmHg) | Leg: Left _____ Right _____   |   |
| <input type="checkbox"/> STRONG (30-40mmHg)  | Foot: Left _____ Right _____  |   |
|  | Knee: Left _____ Right _____<br><small>(only available in Lite and Strong. If Custom, please choose: standard or wraparound.)</small> |   |