



Patient Last Name: _____ Patient First Name: _____
Fitter Last Name: _____ Fitter First Name: _____
Fitter Title: _____ (example: PT/OT/PTA)
Date: _____

FarrowWrap® Thigh High OTS and Custom Order Form

Note: All measurements in cm

Measure medially

Measure posteriorly

*Only necessary if Standard Kneepiece is required.

Measure posteriorly Follow leg contour

A=Floor

Straight distance X

g circumference

f circumference

e circumference

d circumference

c circumference

b² circumference

b¹ circumference

b circumference

a¹ circumference

x straight distance

OTS Components		Thighpiece									
		(Only comes with wrap-around kneepiece)									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(G) Groin		60–70 cm		65–75 cm		70–80 cm		75–90 cm		85–100 cm	
(E) Above Knee		47–52 cm		53–59 cm		60–66 cm		67–73 cm		74–80 cm	
(E-G) Short	FW_-O-TS	18 cm		18 cm		18 cm		18 cm		18 cm	
(E-G) Regular	FW_-O-TR	23 cm		23 cm		23 cm		23 cm		23 cm	
(E-G) Tall	FW_-O-TT	29 cm		29 cm		29 cm		29 cm		29 cm	
		Legpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(C) Mid-Calf		36–43 cm		42–50 cm		48–58 cm		53–63 cm		58–68 cm	
(B) Ankle		21–25 cm		25–30 cm		30–36 cm		36–42 cm		42–50 cm	
(A-D) Regular	FW_-O-LR	35 cm		37 cm		39 cm		41 cm		41 cm	
(B-D)		28 cm		30 cm		31 cm		33 cm		33 cm	
(A-D) Tall	FW_-O-LT	40 cm		42 cm		44 cm		46 cm		46 cm	
(B-D)		33 cm		35 cm		36 cm		38 cm		38 cm	
		Footpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(A¹) Mid-Foot		22–24 cm		25–27 cm		28–30 cm		31–34 cm		35–40 cm	
(X) Regular	FW_-O-FR	16.5 cm		18 cm		19.5 cm		21 cm		23 cm	
(X) Long	FW_-O-FL	18.5 cm		20 cm		21.5 cm		23 cm		25 cm	

Options	
FABRIC	✓
Classic	
LITE	
STRONG	
ThinGrip*	

Use Nonstandard Liner Option (see liner form)	
When a footpiece, legpiece and thighpiece are purchased together, 1 pair of TG® Soft liners are included automatically free. Other liners may cost extra. See liner form & pricing for details. *ThinGrip leg & thighpieces are Custom only. Footpieces available in OTS & Custom.	

CustomComponent	Thighpiece	Wrap-around Standard	Legpiece	FW_-C-L	Footpiece	FW_-C-F	Extend footpiece length to base of toes
Quantity	Right:	Left:	Right:	Left:	Right:	Left:	L R BL: