




Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

## FarrowWrap® Trim-To-Fit Arm Order Form



### FarrowWrap® LITE TTF Arm

| FarrowWrap® TTF Arm                       |            |     | Medium     | XLarge     |
|---|------------|-----|------------|------------|
| (G) axilla circumference                  |            |     | 25 – 44 cm | 30 – 48 cm |
| (F) widest bicep circumference            |            |     | 25 – 42 cm | 25 – 45 cm |
| (E) elbow crease circumference            |            |     | 20 – 38 cm | 25 – 42 cm |
| (D) widest forearm circumference          |            |     | 20 – 36 cm | 20 – 40 cm |
| (C) wrist crease circumference            |            |     | 15 – 31 cm | 15 – 35 cm |
| (E-G) elbow crease to axilla length       |            |     | 15 – 23 cm | 20 – 28 cm |
| (C-E) wrist to elbow crease circumference |            |     | 20 – 25 cm | 20 – 28 cm |
| <b>Left TTF Arm</b>                       | FWLT-T-ATL | Qty |            |            |
| <b>Right TTF Arm</b>                      | FWLT-T-ATR | Qty |            |            |

**This garment only requires 1 hand to don and doff.**

**FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. A free pair of liners is included with each TTF armpiece. Only silver liners are available for this garment.**

**THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY**