




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ (example: PT/OT/PTA)
Date: _____

FarrowWrap® Trim-To-Fit Arm Order Form



FarrowWrap® LITE TTF Arm

FarrowWrap® TTF Arm			Medium	XLarge
(G) axilla circumference			25 – 44 cm	30 – 48 cm
(F) widest bicep circumference			25 – 42 cm	25 – 45 cm
(E) elbow crease circumference			20 – 38 cm	25 – 42 cm
(D) widest forearm circumference			20 – 36 cm	20 – 40 cm
(C) wrist crease circumference			15 – 31 cm	15 – 35 cm
(E-G) elbow crease to axilla length			15 – 23 cm	20 – 28 cm
(C-E) wrist to elbow crease circumference			20 – 25 cm	20 – 28 cm
Left TTF Arm	FWLT-T-ATL	Qty		
Right TTF Arm	FWLT-T-ATR	Qty		

This garment only requires 1 hand to don and doff.

FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. A free pair of liners is included with each TTF armpiece. Only silver liners are available for this garment.

THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY