




Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
Date: \_\_\_\_\_

## FarrowWrap® Off-The-Shelf Arm



FarrowWrap® OTS Arm			Small	Medium	Large
(G) axilla circumference			22 – 31 cm	29 – 39 cm	32 – 45 cm
(E) elbow crease circumference			20 – 27 cm	25 – 34 cm	30 – 40 cm
(C) wrist crease circumference			14 – 18 cm	16 – 21 cm	19 – 25 cm
(C-G) Short	40 – 43 cm	Qty			
(C-G) Regular	44 – 47 cm	Qty			
(C-G) Long	48 – 51 cm	Qty			

**Left or Right must be selected in the “Options” box.**  
**Liners will not be issued with this garment.**

Options	
<b>SIDE</b>	✓
Left	
Right	
<b>Color</b>	✓
Tan	
Black	

*If no color is selected, Tan is the default option.*

### Directions for the Off-The-Shelf Arm

For length, measure from the outside flexor surface of the wrist to the axilla with a slight bend in the arm. Measure circumference at the wrist, elbow crease and axilla. Compare measurements to the size chart to determine correct size. All circumference measurements must correspond to one off-the-shelf size in order for the garment to fit. This garment can be donned and doffed with one hand. **Liners will not be issued with this garment. The OTS Arm is available in LITE fabric only.**

Other: \_\_\_\_\_

**THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY**