




Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

FarrowWrap® Off-The-Shelf Arm



FarrowWrap® OTS Arm		Small	Medium	Large
(G) axilla circumference		22 – 31 cm	29 – 39 cm	32 – 45 cm
(E) elbow crease circumference		20 – 27 cm	25 – 34 cm	30 – 40 cm
(C) wrist crease circumference		14 – 18 cm	16 – 21 cm	19 – 25 cm
(C-G) Short 40 – 43 cm	Qty			
(C-G) Regular 44 – 47 cm	Qty			
(C-G) Long 48 – 51 cm	Qty			

Options	
SIDE	✓
Left	
Right	
Color	✓
Tan	
Black	

**Left or Right must be selected in the “Options” box.
 Liners will not be issued with this garment.**

If no color is selected, Tan is the default option.

Directions for the Off-The-Shelf Arm

For length, measure from the outside flexor surface of the wrist to the axilla with a slight bend in the arm. Measure circumference at the wrist, elbow crease and axilla. Compare measurements to the size chart to determine correct size. All circumference measurements must correspond to one off-the-shelf size in order for the garment to fit. This garment can be donned and doffed with one hand. **Liners will not be issued with this garment. The OTS Arm is available in LITE fabric only.**

Other: _____

THIS IS FOR THE INTENDED USE OF LUNA MEDICAL ONLY