



Bellavar® and Custom Seamless Soft Order Form

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

Product / Brand	Quantity		Sand	Black	Bronze	Caramel	Navy	Cranberry	Espresso	Sun Bronze
	Left	Right								
Seamless Soft 18-21 mmHg* (CCL 1)										
Seamless Soft 23-32 mmHg* (CCL 2)										
Seamless Soft 34-46 mmHg* (CCL 3)										
Bellavar® 23-32 mmHg* (CCL 2)										
Bellavar® 34-46 mmHg* (CCL 3)										

Basic Styles:

AD AF AG AG-T AG-HT AT

Options:

Closed toe Open toe Short foot (closed)


Special Options:

AD	<input type="checkbox"/> No Silicone	<input type="checkbox"/> Silicone dotted band 2.5 cm
	<input type="checkbox"/> Silicone dotted band 5 cm	<input type="checkbox"/> SoftFit™ (only in CCL1 & CCL2)***

AF/AG	<input type="checkbox"/> No Silicone	<input type="checkbox"/> Silicone dotted band 5 cm
	<input type="checkbox"/> Silicone lace band 6 cm	<input type="checkbox"/> Silicone Soft band 6 cm**
	<input type="checkbox"/> Sensitive Band (Seamless Soft Only)	

AT	<input type="checkbox"/> Maternity	<input type="checkbox"/> Fly for Men
	<input type="checkbox"/> Full compression	<input type="checkbox"/> Regular Adjustable Waist band
	<input type="checkbox"/> Waist band 2.5 cm**	<input type="checkbox"/> Waist band 5.0 cm**
	<input type="checkbox"/> Open Pubis	<input type="checkbox"/> Mesh Crotch

Form 57021 must accompany this form.			
Circum.(c)	Length (l)	Length (l)	
cT	K2-T	lT	
cH	K1-T	lH	
Circumference (c)		Length (l): Taken from each landmark to floor.	
Left	Right	Left	Right
cG		lG	
cF		lF	
cE		lE	
cD		lD	
cC		lC	
cB1		lB1	
cB		lB	
cY		lZ (closed toe)	
cA		lA (open toe)	



Foot length open toe lA _____ Foot length closed toe lZ _____

(Not available in slant open or slant closed toe, only straight.)

Comments: _____

*Design Pressure **Not available in Full Compression or Bellavar* ***Not available in Bellavar*

.....Take measurements on edema-free extremities only. All measurements must be recorded in cm.