



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

JOBST NON-CUSTOM HAND GLOVE/GAUNTLET MEASUREMENT FORM AND SIZING CHART



Measure in centimeters

Left ☐ Right ☐

Measurement Points:

Circumference **point 'a'** is found across the metacarpus.

Circumference **point 'c'** is found distally on the arm where the wrist meets the hand.

One should find the smallest circumference at this point.

Circumference **'a'** = _____

Circumference **'c'** = _____

Style	
Gauntlet w/ thumb stub only	☐
Glove	☐

Material	Class 1 15-20 mmHg	Class 2 20-30 mmHg	Color Available
Bella Strong Glove	☐	☐	Beige ☐ Black ☐
Bella Strong Gauntlet	☐	☐	Beige ☐ Black ☐
Bella Lite Glove	☐	☐	Beige ☐ Black ☐
Bella Lite Gauntlet	☐	☐	Beige ☐ Brown ☐ Black ☐ *Brown and Black only available in 20-30 mmHg

Sizing chart:

JOBST® Bella™ Strong Gauntlet				
	WRIST		PALM	
SIZE	MIN.	MAX.	MIN.	MAX.
1	5 1/2" (14 cm)	7" (18 cm)	6" (15 cm)	6 3/4" (17 cm)
2	5 1/2" (14 cm)	7" (18 cm)	6 3/4" (17 cm)	7 1/2" (19 cm)
3	7" (18 cm)	8 1/4" (21 cm)	6 3/4" (17 cm)	7 1/2" (19 cm)
4	5 1/2" (14 cm)	7" (18 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)
5	7" (18 cm)	8 1/4" (21 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)
6	6 1/4" (16 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)	9" (23 cm)
7	7 1/2" (19 cm)	8 1/2" (22 cm)	8 1/4" (21 cm)	9" (23 cm)
8	7 1/2" (19 cm)	8 1/4" (22 cm)	9" (23 cm)	10 1/4" (26 cm)

Bella Strong Glove

Size	cA	cC
1	16.6 cm - 18.0 cm	15.6 cm - 16.7 cm
2	16.6 cm - 18.0 cm	16.7 cm - 18.1 cm
3	18.0 cm - 19.5 cm	15.9 cm - 17.0 cm
4	18.0 cm - 19.5 cm	17.0 cm - 18.5 cm
5	18.0 cm - 19.5 cm	18.5 cm - 20.4 cm
6	19.5 cm - 21.5 cm	17.2 cm - 18.7 cm
7	19.5 cm - 21.5 cm	18.7 cm - 20.6 cm
8	21.5 cm - 23.1 cm	18.3 cm - 20.1 cm


JOBST® Bella™ Lite Gauntlet		
	PALM	WRIST
S	5 7/8" - 7 3/8" (15-19 cm)	5 3/8" - 7 5/8" (14-20 cm)
M	7 1/2" - 8 5/8" (19-22 cm)	6 1/8" - 8 1/2" (16-22 cm)
L	8 3/4" - 10 1/4" (22-26 cm)	6 7/8" - 9 1/4" (18-24 cm)

Bella Lite Glove

Measurements Table

Circumference in cm:

Size	1	2	3	4	5	6	7	8
A	16.6-18.0 cm	16.6-18.0 cm	18.0-19.5 cm	18.0-19.5 cm	18.0-19.5 cm	19.5-21.5 cm	19.5-21.5 cm	21.5-23.1 cm
B	15.6-16.7 cm	16.7-18.1 cm	15.9-17.0 cm	17.0-18.5 cm	18.5-20.4 cm	17.2-18.7 cm	18.7-20.6 cm	18.3-20.1 cm



Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1057 W. Grand Ave · Suite 1 · Chicago, IL 60642 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com