




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



Classic Arm Sleeve AG1 Ready-to-Wear

SHIPPING ADDRESS			
Patient Name			
Address			
Attention			
City		State	
Phone		Zip	



Polartec [®] Power Dry [®] Colors & Channeling	Side	Size		Length
<input type="checkbox"/> Black - Slimline Channeling	<input type="checkbox"/> Right Arm	<input type="checkbox"/> XSmall	<input type="checkbox"/> XSmall MAX	<input type="checkbox"/> XShort
<input type="checkbox"/> Royal Blue - Standard Channeling	<input type="checkbox"/> Left Arm	<input type="checkbox"/> Small	<input type="checkbox"/> Small MAX	<input type="checkbox"/> Short
<input type="checkbox"/> Royal Blue - Slimline Channeling		<input type="checkbox"/> Medium	<input type="checkbox"/> Medium MAX	<input type="checkbox"/> Regular
<input type="checkbox"/> JoViJacket - Black Nylon & Spandex Powernet (JoViJackets are recommended as they provide the additional compression needed for maximum fit & effectiveness.)		<input type="checkbox"/> Large	<input type="checkbox"/> Large MAX	<input type="checkbox"/> Long
		<input type="checkbox"/> XLarge	<input type="checkbox"/> XLarge MAX	<input type="checkbox"/> XLong
Comments/Quantity:				
Additional Items				
<input type="checkbox"/> Dycem [®] - donning aid				
<input type="checkbox"/> Easy Slide - donning aid				
<input type="checkbox"/> Prepaid Reduction Option				



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____



Classic Arm Sleeve AG1 Ready-to-Wear

All measurements are provided in centimeters.

Classic Arm Sleeve AG1 Size Chart			
	Least Wrist Circumference	Least Elbow Circumference	Axilla Circumference
	"C" Landmark	"E" Landmark	"G" Landmark
XSmall	14 - 15.5	20 - 23	23 - 27.5
XSmall MAX		22 - 25	27 - 31
Small	15.5 - 16.5	23 - 26	27.5 - 32
Small MAX		25 - 28	31 - 36
Medium	16.5 - 18	26 - 29	32 - 36
Medium MAX		28 - 31	36 - 40
Large	18 - 19	29 - 33	36 - 41
Large MAX		31 - 35	40 - 44.5
XLarge	19 - 20.5	33 - 37.5	41 - 45
XLarge MAX		35 - 39.5	44.5 - 49

MAX sizes are 2cm larger at the elbow and 4cm larger at the axilla.

Classic Arm Sleeve AG1 Length Chart	
Wrist to Axilla (C to G) (measured with arm hanging down at side)	
XShort	38
Short	40.5
Regular	43
Long	45.5
XLong	48