




Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____



Classic Lower Leg AD Ready-to-Wear



Polartec® Power Dry® Colors	
<input type="checkbox"/>	Black
<input type="checkbox"/>	Royal Blue

<input checked="" type="checkbox"/>	JoViJacket - Black Nylon & Spandex Powernet (JoViJackets are recommended as they provide the additional compression needed for maximum fit & effectiveness.)
-------------------------------------	---

Comments/Quantity:

Size	Length
<input type="checkbox"/> Small	<input type="checkbox"/> Regular
<input type="checkbox"/> Medium	
<input type="checkbox"/> Large	
<input type="checkbox"/> XLarge	
<input type="checkbox"/> 2XLarge	

Additional Items
<input type="checkbox"/> Dycem® - donning aid
<input type="checkbox"/> Arion Easy Slide - donning aid

Dycem® is a registered trademark of Dycem Ltd.



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



Classic Lower Leg AD Ready-to-Wear

All measurements are in centimeters (cm).

Classic Lower Leg AD Size Chart			
	Least Ankle Circumference	Widest Calf Circumference	Least Knee Circumference
	"B" Landmark	"C" Landmark	"D" Landmark
Small	20.5-22	35.5-38	32-34.5
Medium	22-23	38-40.5	34.5-37
Large	23-24	40.5-43	37-39.5
XLarge	24-25.5	43-45.5	39.5-42
2XLarge	25.5-27	45.5-48	42-44.5

Classic Lower Leg AD Length Chart	
Floor to Tibial Tuberosity (A to D) (measured medially)	
Regular	42



Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Shoe Size	BNR	UOM / Box	Order Qty.
X-Small	≤2	7965803	1	
Small	2.5-5.5	7965804	1	
Medium	6-8	7965802	1	
Large	8.5-11	7965902	1	
X-Large	≥11.5	7966001	1	