



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

## JUZO NON-CUSTOM HAND GLOVE/GAUNTLET MEASUREMENT FORM AND SIZING CHART

**Measure in centimeters**

Left  Right

**Measurement Points:**

Circumference point 'a' is found across the metacarpus.  
 Circumference point 'c' is found distally on the arm where the wrist meets the hand.  
 One should find the smallest circumference at this point.

Circumference 'a' = \_\_\_\_\_

Circumference 'c' = \_\_\_\_\_

Style	
Gauntlet w/ thumb stub only	<input type="checkbox"/>
Glove w/ finger stubs	<input type="checkbox"/>


Material	Class 1 18-21 mmHg	Class 2 23-32 mmHg	Color Available
Expert <small>*formerly Helastic</small>	3021 <input type="checkbox"/>	3022 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/>

Material	Class I 20-30 mmHg	Class II 30-40 mmHg	Colors Available
2300 Series Gauntlet	2301 <input type="checkbox"/>	2302 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/> Finger Stubs <input type="checkbox"/>
Dreamsleeve Soft 2000 Gauntlet	2001 <input type="checkbox"/>	2002 <input type="checkbox"/>	Beige <input type="checkbox"/> Black <input type="checkbox"/> Chestnut <input type="checkbox"/> Cinnamon <input type="checkbox"/> <small>_____ <input type="checkbox"/> * Seasonal colors available</small>

\*Soft 2000 are only available as gauntlet with thumb stub.

\* To view seasonal colors, visit [www.dreamsleeve.com](http://www.dreamsleeve.com)

### Size Chart for Juzo® 3021/3022 Gauntlet



	1		2		3		4		5		6	
	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm
①	6.5	5	6.5	5	6.5	5	7.5	6	7.5	6	7.5	6
②	6.2	5	6.7	5.5	7.1	5.9	7.6	6.4	8	6.8	8.5	7.3
③	6	5	6.3	5.3	6.7	5.7	7.1	6	7.5	6.4	7.6	6.7
④	5.7	4	6	4.2	6.3	4.5	6.7	4.8	7.1	5.1	7.4	5.4
⑤	5.5	4	5.9	4.2	6.3	4.5	6.6	4.8	7	5.1	7.3	5.4
a	17.5 cm 6 7/8"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"		23.5 cm 9 1/4"		25 cm 9 3/4"	
c	14.5 cm 5 3/4"		16 cm 6 1/4"		17.5 cm 7"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"	

### Size Chart for Juzo® 1102, 2001 / 2002 Gauntlet

	1-S	2-M	3-L	4-XL
a	17.5-19.5 cm 6 7/8"-7 5/8"	19.5-21.0 cm 7 5/8"-8 1/4"	21.0-23.0 cm 8 1/4"-9 1/8"	23.0-25.5 cm 9 1/8"-10 1/8"
c	14.5-16.5 cm 5 5/8"-6 1/2"	16.5-18.0 cm 6 1/2"-7"	18.0-20.0 cm 7"-7 7/8"	20.0-22 cm 7 7/8"-8 5/8"

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1057 W. Grand Ave · Suite 1 · Chicago, IL 60642 · Phone (800) 380-4339 · Fax (888) 696-0299 · [www.lunamedical.com](http://www.lunamedical.com)