



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

## 20-30 and 30-40 mmHg Compression Stockings 2501/02, 2581/82, 2001/02, 2081/82, 4601/02 3511/12, 2061/62

	I	II	III	IV	V
<b>cG</b>	41-60 cm 16 1/4"-23 1/2"	50-68 cm 19 1/2"-26 3/4"	54-75 cm 21 1/4"-29 1/2"	57-79 cm 22 1/2"-31"	62-85 cm 24 1/2"-33 1/2"
<b>cG</b> with silicone	45-55 cm 17 3/4"-21 1/2"	50-61 cm 19 1/2"-24"	54-68 cm 21 1/4"-26 3/4"	57-73 cm 22 1/2"-28 3/4"	62-78 cm 24 1/2"-30 3/4"
<b>cC</b>	29-38 cm 11 1/2"-15"	34-43 cm 13 1/4"-17"	37-49 cm 14 1/2"-19 1/4"	41-53 cm 16 1/4"-21"	46-58 cm 18"-23"
<b>cB</b>	18-21 cm 7"-8 1/4"	21-24 cm 8 1/4"-9 1/2"	24-27 cm 9 1/2"-10 3/4"	27-31 cm 10 3/4"-12 1/4"	31-35 cm 12 1/4"-13 3/4"

### length A-D

**A-D Regular** - 15 3/4-18 1/4" / 40-46 cm

**A-D Short** - 13-15 3/4" / 34-40 cm

**A-D Petite** - < 13" / 33 cm

### length A-G

**A-G / A-T Regular** - 28 1/4-32 3/4" / 72-83 cm

**A-G / A-T Short** - 24 3/4-28 1/4" / 63-72 cm

**A-G / A-T Petite** - 21 3/4-24 3/4" / 55-63 cm

PCSZ-01-02c

## Size Chart for Juzo® 40-50 mmHg Compression Stockings and Pantyhose

	I	II	III	IV	V	VI
<b>cG</b>	43-56 cm 17"-22"	47-62 cm 18 1/2"-24 1/2"	49-66 cm 19 1/4"-26"	55-73 cm 21 3/4"-28 3/4"	59-77 cm 23 1/4"-30 1/4"	64-82 cm 25 1/4"-32 1/4"
<b>cG</b> with silicone	46-53 cm 18"-21"	51-58 cm 20"-23"	53-62 cm 21"-24 1/2"	60-69 cm 23 1/2"-27 1/4"	64-73 cm 25 1/4"-28 3/4"	69-78 cm 27 1/4"-30 3/4"
<b>cC</b>	27-34 cm 10 1/2"-13 1/2"	29-37 cm 11 1/2"-14 1/2"	33-41 cm 13"-16 1/4"	36-45 cm 14 1/4"-17 3/4"	39-47 cm 15 1/2"-18 1/2"	42-50 cm 16 1/2"-19 3/4"
<b>cB</b>	18-20 cm 7"-7 3/4"	20-22 cm 7 3/4"-8 1/2"	22-24.5 cm 8 1/2"-9 1/2"	24.5-27 cm 9 1/2"-10 1/2"	27-30 cm 10 1/2"-11 3/4"	30-33 cm 11 3/4"-13"

### length A-D

**A-D Regular** - 15 3/4-18 1/4" / 40-46 cm

**A-D Short** - 13-15 3/4" / 34-40 cm


**A-D Petite** - < 13" / 33 cm

### length A-G

**A-G Regular and A-T Regular** - 28 1/4-32 3/4" / 72-83 cm

**A-G Short and A-T Short** - 24 3/4-28 1/4" / 63-72 cm

**A-G Petite and A-T Petite** - 21 3/4-24 3/4" / 55-63 cm



**IA**

**hip circumference**

**A-T Form A** < 41" / 104 cm

**A-T Form B** > 41" / 104 cm