




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

Size Charts For Dynamic Max Compression Stockings 20-30 mmHg and 30-40 mmHg



20-30, 30-40 mmHg Dynamic Max Compression Stockings 3511 / 3512

	III	IV	V	
cG	60-82 cm 23 3/4"-32 1/4"	62-87 cm 24 1/2"-34 1/4"	77-94 cm 30 1/4"-37"	length A-D
				A-D Regular - 15 3/4"-18 1/4" / 40-46 cm A-D Short - 13-15 3/4" / 34-40 cm
cC	49-55 cm 19 1/4"-21 3/4"	53-60 cm 21 1/4"-23 1/2"	55-66 cm 21 3/4"-26"	length A-G
				A-G Regular - 28 1/4"-32 3/4" / 72-83 cm A-G Short - 24 3/4"-28 1/4" / 63-72 cm
cB	24-27 cm 9 1/2"-10 3/4"	27-32 cm 10 3/4"-12 3/4"	32-37 cm 12 3/4"-14 3/4"	

PCSZ-01-47