



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

Instructions for precise measuring of Juzo® thorax compression segments

Circumferences of thorax:

- cT** = Waist measurement
- cN** = Chest measurement under the armpit
- cH** = Arm circumference, measured from the shoulder straight down under the armpit
- cK** = Hip circumference, in the case of the extended thorax bandage at least as low as the hip bones

Length measurements of thorax (to be measured straight down):

To be recorded on the measurement schematic as per the arrow lines.

- l TS** = from the shoulder to the waist
- l KT** = from the waist to the hips
- l KK'** = from the height of the hips at the front down through the crotch to the height of the hips at the back
- l QU** = width of the neck cutout at the front
- l RS** = depth of the neck cutout at the front

If you do not want the back part closed all the way to the back of the neck, please specify the depth of the rear cutout. The width of the cutout at the back will be the same as the width of the neck cutout.

If **chest cutouts or knitted-in cups** are required, please specify the cup size.

No additional inelastic bra is required for compression of the chest.

The knitted-in cups are knitted from a virtually inelastic thread and therefore have a slightly compressive effect.

The thorax bandage can be connected to a compression stocking by means of a hook and loop fastener.

Where this is required, send both measurement schematics and include a suitable remark.

The **inelastic back zone** is an approx. 3 cm wide area of the back that has been knitted with an inelastic fiber.

It avoids the adoption of a relieving posture in children.

In the style with **external seams**, all seams are located on the outside to ensure an even flatter seam on the inside.

To achieve an optimum pressure in the chest area, we recommend using arm extensions (T-shirt style) or compression sleeves.

Instructions for precise measuring of associated compression sleeves or arm extensions

Circumferences:

- cC** = Wrist
- cC'** = 6 cm above "cC" for adults
individually selected for children
- cD** = Mid forearm
- cE** = Mid elbow
- cF** = Mid upper arm
- cG** = Height of the armpit

Length measurements:

All length measurements are taken from "C" or, in the case of arm extensions, from the relevant end measurement.

l GH is to be measured from the height "cG" to the shoulder.



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

Quantity: piece(s)

Please check: <input checked="" type="checkbox"/>	Compression classes	
	1	2
with seam (flat knitted)		
Juzo® Expert Colors: <input type="checkbox"/> Beige <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Grey	<input type="checkbox"/> 3021	<input type="checkbox"/> 3022
<input type="checkbox"/> Dark blue <input type="checkbox"/> Mocca <input type="checkbox"/> Black		
Juzo® Expert Silver (Color: Beige)	<input type="checkbox"/> 3021	<input type="checkbox"/> 3022

Area of application:

- Lymphology Scar treatment

Model:

- Opening: Mid front Mid back
- With zipper With hook & eye fastener With velcro fastener

- Slip on
- With arm sleeve
- Without arm sleeve
- Breast opening, cup size
- Breast cup seamless, cup size

- Stand up collar
- Neck circumference: cm

Stand up collar height: cm
(Measurements IQU and IRS are not needed)

- Attached on a body part of a compression pantyhose (for velcro closure at "T" please enclose the corresponding order form)
- Adhesive border at "T"
- Inelastic back zone
- External seams
- With gusset (length of gusset: 28 cm, width: 10 cm)
- Body fastener (gusset with hook fastener, only in conjunction with fastening option) · Body height: cm
- Anatomically flexed middle of elbow

Special requests:

.....

.....

.....

.....

From:

.....

.....

.....

.....

